

# New at Wolf X-Ray

## Image Plate Cleaning Wipes



Clean Image makes taking care of imaging plates quick, easy and affordable. Treated with a proprietary Methyl alcohol formula, Clean Image wipes easily remove fingerprints, dirt and artifacts. Simply tear open a packet and wipe down your plates. Wipes are individually wrapped and disposable, making them perfect for single-use.

Clean Image Wipes **#58150**  
Packaged 50/box

**Wolf**<sup>™</sup> WOLF X-RAY CORPORATION  
100 West Industry Court, Deer Park, NY 11729  
1-800-356-9729 • 631-242-9729 • Fax 631-242-1001  
X-RAY CORPORATION [www.WolfXray.com](http://www.WolfXray.com) • [info@WolfXray.com](mailto:info@WolfXray.com)