

Wolf™

PROTECTIVE APPAREL AND RELATED EQUIPMENT



Easy Wrap Apron

Our Easy Wrap Apron is the professional standard. The two side closures allow you to custom fit the shoulder weight distribution. A handy pocket is sewn on the breast.

.5 mm LE Protection.

Size		Standard Lead	Light Weight Lead	LEAD FREE
Small	22 x 37"	65025	65025LW	65025TB
Medium	24 x 40"	65021	65021LW	65021TB
Large	24 x 42"	65023	65023LW	65023TB
Extra Large	25 x 40"	65026	65026LW	65026TB
Double Extra Large	26 x 42"	65027	65027LW	65027TB

With Collar

Size		Light Weight Lead	LEAD FREE
Small	22 x 37"	65025TC	65025TC-TB
Medium	24 x 40"	65021TC	65021TC-TB
Large	24 x 42"	65023TC	65023TC-TB
Extra Large	25 x 40"	65026TC	65026TC-TB
Double Extra Large	26 x 42"	65027TC	65027TC-TB



Easy Wrap with Collar

Attached to the apron with strap for easy removal.

Coat Apron

The Coat Apron features a 2 inch wide adjustable belt that helps "lift" the apron weight off your shoulders. The open back design allows you great freedom of movement.

.5 mm LE Protection.

Size		Standard Lead	Light Weight Lead	LEAD FREE
Small	22 x 37"	63073	63073LW	63073TB
Medium	24 x 40"	63001	63001LW	63001TB
Large	24 x 42"	63011	63011LW	63011TB
Extra Large	25 x 40"	63012	63012LW	63012TB
Double Extra Large	26 x 42"	63013	63013LW	63013TB

With Collar

Size		Light Weight Lead	LEAD FREE
Small	22 x 37"	63073TC	63073TC-TB
Medium	24 x 40"	63001TC	63001TC-TB
Large	24 x 42"	63011TC	63011TC-TB
Extra Large	25 x 40"	63012TC	63012TC-TB
Double Extra Large	26 x 42"	63013TC	63013TC-TB

Place the Color/Pattern # next to the Apron stock # when ordering



Coat Apron with Collar

attached to the apron with a strap for easy removal.



Apron and Vest

A unique protective apron design that provides complete front and back coverage with the apron weight distributed evenly across the shoulders and hips. The apron is a wrap around style which ensures a perfect fit.

Now Wolf's popular Apron and Vest is available specifically sized for Women and Men (see stock number tables below) and is also available with an optional Thyroid Collar. The detachable Thyroid Collar is connected to the vest by a strap for easy removal.

LE Protection, front .5 mm and back .3 mm.

Women's Style Stock Numbers

Size	Apron: L x W	Vest: L x W	Standard Lead	Light Weight Lead	LEAD FREE
Small	19 x 30-34	21 x 19	66077	66077LW	66077TB
Medium	20 x 35-40	22 x 21	66078	66078LW	66078TB
Large	21 x 41-46	23 x 23	66079	66079LW	66079TB
Extra Large	22.5 x 47-50	24 x 25	66080	66080LW	66080TB

Women's Style With Collar

Size	Apron L x W	Vest L x W	Light Weight Lead	LEAD FREE
Small	19 x 30-34	21 x 19	66077TC	66077TC-TB
Medium	20 x 35-40	22 x 21	66078TC	66078TC-TB
Large	21 x 41-46	23 x 23	66079TC	66079TC-TB
Extra Large	22.5 x 47-50	24 x 25	66080TC	66080TC-TB

Men's Style Stock Numbers

Size	Apron: L x W	Vest: L x W	Standard Lead	Light Weight Lead	LEAD FREE
Small	20 x 35-40	24 x 21	66081	66081LW	66081TB
Medium	21 x 41-46	25 x 23	66082	66082LW	66082TB
Large	22.5 x 47-50	26 x 25	66083	66083LW	66083TB
Extra Large	23 x 51-56	27 x 27	66084	66084LW	66084TB

Men's Style With Collar

Size	Apron: L x W	Vest: L x W	Light Weight Lead	LEAD FREE
Small	20 x 35-40	24 x 21	66081TC	66081TC-TB
Medium	21 x 41-46	25 x 23	66082TC	66082TC-TB
Large	22.5 x 47-50	26 x 25	66083TC	66083TC-TB
Extra Large	23 x 51-56	27 x 27	66084TC	66084TC-TB

Apronette



Wolf's classic protective apron is held on by a belt with a velcro closure.

.5 mm LE Protection.

Size	Stock #
12 x 12"	69095
12 x 18"	69096
12 x 24"	69097
18 x 24"	69098



Lumbar Support Belt

The perfect accessory for any Wolf Apron, the Lumbar Support Belt custom fits like a glove and provides hours of comfortable support. Simply wrap the belt around your waist with the velcro closure. Pull the tab and feel the cables in the back of the belt custom hug your lumbar region. When you've reached your individual comfort level, secure the tab to the belt.




Size	Part No.
Small	61005
Medium	61006
Large	61007
X Large	61008

O. R. Quick Drop Apron

Expressly for maintaining an uncontaminated sterile field. Made to be worn between a scrub-suit and sterile gown, all the surgeon has to do is flip back the two shoulder tabs and the apron falls to the ground. Thyroid Collar optional.
.5 mm LE Protection.

Size		Light Weight Lead	LEAD FREE
Small	22 x 37"	68094LW	68094TB
Medium	24 x 38"	68098LW	68098TB
Large	24 x 42"	68099LW	68099TB
Extra Large	25 x 40"	68100LW	68100TB
Double Extra Large	26 x 42"	68101LW	68101TB



Special Procedure Apron

The Special Procedure has an elastic belt that helps lift the apron weight off your shoulders. A full back panel gives you all around protection and the over-the-shoulder velcro closure makes putting it on and removing it a cinch.


LE Protection, front .5 mm and back .3 mm.

Size		Standard Lead	Light Weight Lead	LEAD FREE
Small	21 x 36"	67086	67086LW	67086TB
Medium	22 x 38"	67087	67087LW	67087TB
Large	24 x 38"	67088	67088LW	67088TB
Extra Large	24 x 40"	67089	67089LW	67089TB
Double Extra Large	25 x 42"	67090	67090LW	67090TB



Gonad Shield

.5 mm LE Protection
 Has a belt. Beige only
 #75500-23



3 Gonad Shield Set

.5 mm LE Protection
 S 3 x 5", M 5 x 7", L 7 x 9".
 Navy blue only.
 #75501-22



SofTex Thyroid Collar

A soft pliable collar with velcro closure,
.5 mm LE Protection
 Comes in all colors.
 # 75070 + color #



Thyroid Collar

Held in place by velcro, comes in all colors.
.5 mm LE Protection
 #75060 + color #



Front Closing Special Procedure Apron

The Front Closing Special Procedure Apron allows weight to be distributed across the shoulders and comfort is enhanced by an adjustable belt. The Front Closing Special Procedure Apron has a full back panel for all-around protection. The apron is held closed by three velcro tabs. *LE Protection, front .5 mm and back .3 mm.*

Size		Standard Lead	Light Weight Lead	LEAD FREE
Small	21 x 36"	68091	68091LW	68091TB
Medium	22 x 38"	68093	68093LW	68093TB
Large	24 x 38"	68095	68095LW	68095TB
Extra Large	24 x 40"	68096	68096LW	68096TB
Double Extra Large	25 x 42"	68097	68097LW	68097TB

The Conventional Apron

Simple, light frontal protection, our Conventional Apron is the professional standard. The apron is held in place by two adjustable straps that cross in the back to assure a snug fit. This apron is ideal as a veterinary apron. *.5 mm LE Protection.*




Size		Standard Lead	Light Weight Lead	LEAD FREE
Small	22 x 33"	62008	62008LW	62008TB
Medium	24 x 39"	62018	62018LW	62018TB
Large	24 x 42"	62019	62019LW	62019TB
Extra Large	25 x 40"	62020	62020LW	62020TB
Double Extra Large	26 x 42"	62021	62021LW	62021TB

Wolf Apron Color / Pattern Choices




11 Rose Pink 12 Purple 14 Teal 15 Powder Blue 16 Medium Blue 18 Burgundy 19 Grey



20 Royal Blue 22 Navy Blue 23 Beige 25 Black 35 Red 36 Hunter Green



30 Tye Dye



32 Hawaiian



33 Camouflage



34 Dog Paws

Place the Color/Pattern # next to the Apron stock # when ordering

Accessories



Superflex Protective Gloves

Molded lead glove covered in Naugahyde.

Length		Equiv.	Pair	R Only	L Only	Cotton Inserts
12 inch	Reg.	.25mm	12420	12520	12620	
12 inch	Reg.	.5mm	12419	12519	12619	12001-02
15 inch	Long	.5mm	12422	12522	12622	
15 inch	Long	.25mm	12426	12526	12626	
15 inch	Flared	.5mm	12425	12525	12625	



Comfort Flex Gloves

Perfect for vet. use. Small 12406-S, Medium 12406-M, Large 12406-L

Sensi Flex Gloves

For Cath/EP labs, Ortho, Endoscopy, etc. Available sizes 7 to 8 - 1/2



Attenuation

60 KVP	80 KVP	100 KVP	120 KVP	Lead Equiv.	Stock No.
45%	35%	26%	23%	.22 m	12407
55%	45%	35%	31%	.30 m	12408

Slitted Mittens

A perfect combination of protection and dexterity, our slitted mittens allow you the option of exposing your fingers while protecting your hand.

.5 mm LE Protection. #12410



Deca-Rak

The Deca-Rak has 10 pivoting arms that fan out to 180°. It is mounted on locking casters and has a white finish throughout.

Specifications:
55"H x 51"W x 20"D with center arm at 90° angle from racks depth is 26". Approx. Shipping Weight: 125 lbs.
#16411



Penta-Rak

Wolf's wall mounted 5 apron storage unit has 5 pivoting apron hangers that fan out 180° and has a white finish throughout. The Penta-Rak holds up to 110 lbs. maximum weight.

Specifications: 4 5/8"H x 17"W with center arm at 90° angle from racks depth is 22". Approx. Shipping Weight: 30 lbs.
#16410

Free Standing Penta-Rak

The Penta-Rak is available as a free standing unit. It looks the same as the Deca-Rak (shown left) but is half the length.
#16405



Tri-Rak

Our wall mounted Tri-Rak holds a protective apron and a pair of protective gloves.

Specifications:
16"H x 28"W x 8"D. Approx. Shipping Weight: 9 lbs.
#16400



Uni-Hanger

Hangs any protective apron on a wall hook or in a conventional closet.
#16401



Glove-Rak Kit

With Wolf's Glove-Rak Kit you can add storage for 2 pairs of gloves to a Penta-Rak or 4 pairs of gloves to a Deca-Rak. Glove posts are simply screwed onto the hanger unit via pre-drilled holes.
Note: Only Raks manufactured since 10/04 have pre-drilled holes.
#16415



Skirt Hanger

Our new skirt hanger can be attached to a Deca or Penta-Rak arm as a perfect way to store a protective skirt or Apronette.
#16402



Wolf X-Ray Corporation

100 West Industry Court, Deer Park, NY 11729

Phone: 631-242-9729 • Toll free: 800-356-9729 • Fax: 631-242-1001

email: info@wolfxray.com • www.wolfxray.com