



**INSTRUCTIONS FOR USE
21695, COOL-LITE**

REF

**19103, 19204: ANGULAR SAFELIGHT
19198: SWIVEL LIGHT**



1) Indication for use:

- a. COOL -LITE: This unique viewing device provides concentrated light like our Intensi-Spot, but eliminates the heat by the use of a fluorescent bulb. Energy efficient, the Cool-lite can be left on all day. The Cool-Lite tilts to any angle, has a baked enamel finish, and includes holes in the base for easy wall mounting. The Cool-Lite comes with a 6' power cord with a 3-prong plug.
- b. ANGULAR LIGHT: ideal where a direct light over the developing tank or loading bench is wanted. Easy to install via 2 keyholes, the Angular Safelights GBX Type filter measures 5" x 7". The Angular comes complete with a bulb, and a 6' cord with a three-prong plug. Optional SENSOR turns the Safelight on automatically when regular light is shut off. (The unit can still be shut off permanently with the on/off switch).
- c. SWIVEL LIGHT: Perfect for direct spotlighting this wall mounted unit swivels to any position. Easy to install via 2 keyholes, the Swivel Safelight includes a GBX Type filter, light bulb, and a 6' cord with a three-prong plug. Optional SENSOR turns the Safelight on automatically when regular light is shut off. (The unit can still be shut off permanently with the on/off switch).

COOL-LITE

2) Instructions for Use

Note: It is always recommended to

- A. Tilts to any angle
- B. Note: Cool-Lite can be left on all day.
- C. Holes in the base for wall mounting
- D. Fluorescent Light Replacement Bulb (75 w)



- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) Precautions/Warnings: California's Proposition 65
- 5) Disposal: Product should be disposed of according to applicable local laws.

rev 01, 21695, 19103, 19104, 19198
08/2021