



Functionality	Performance Monitor		VitalsRx Monitor*		O2 Vibe		
MSRP	\$299		\$499		\$199		
20-Second ECG/EKG waveform	✓		✓				
Continuous ECG/EKG monitoring via USB electrodes or patch			✓				
20-Second SpO2 capture	✓		✓		✓		
Continuous SpO2 monitoring via USB finger sensor or cable free			✓		✓		
Overnight sleep screening report			✓		✓		
Heart rate monitoring	✓		✓		✓		
Heart rate variability (HRV)	✓		✓				
Respiration rate (RR)	✓						
Diastolic & systolic blood pressure ¹	BP Cuff	✓	Optional	BP Cuff	✓	Optional	N/A
20-Second systolic blood pressure	✓					N/A	
3-Second thermometer	✓		✓			N/A	
Relax-Me stress relief exercises	✓		✓			N/A	
Pedometer: steps, calories, distance, etc.	✓		✓			Steps	
Medication reminders	✓		✓			N/A	
Profiles or users	8		8			1	
Free iOS & Android mobile application	✓		✓			✓	
FDA Exempt – General Wellness Guidelines July 2016	✓		N/A			N/A	
FDA Cleared – 510(k) No.150869 or Registered	N/A		✓			✓	
CE Mark Certified – CE No. HD 60095648 0001	✓		N/A			✓	

¹ Optional Integrated wrist BP Cuff available for purchase \$99

* **VitalsRx Monitor** For Sale to Healthcare Professionals Only - Rx Required for Consumer Use

All devices have settings for English and Spanish languages.