

MobilDrTech / Riester

# RCS-100



## PRODUCT COMPARISON



Riester  
RCS-100



Horus DCS-100

<b>Image Resolution</b>	3264 x 2448	1920 x 1080
<b>MegaPixels</b>	8MP	2 MP
<b>Operating Time</b>	3.5 hours	3 hours
<b>Charging Time</b>	3 hours	5 hours
<b>File Storage</b>	16GB to 32GB	Supports 2GB to 32GB
<b>Warranty</b>	2 years	1 year
<b>Focus</b>	Manual/Auto Focus	Manual Focus
<b>Display on device</b>	5" full HD colour TFT-LCD	3.5" Full Color TFT LCD
<b>Image Format</b>	JPEG (photo)/ H.264 (Video)	JPEG (photo)/ H.264 (Video)