



# O2 Vibe Monitor

Settings for English and Spanish languages

The Only Sleep Monitor That Vibrates You To Better Sleep



**O2 Vibe** (17+)  
BodiMetrics, LLC  
★★★★★ 2.8, 11 Ratings  
Free



## Features and Metrics

- 20-Second SpO2 Capture
- Continuous SpO2 Monitoring
- Overnight Sleep Screening Report
- Heart Rate
- Pedometer Steps
- FDA Registered / CE Certified

- SpO2/Oxygen
- Heart Rate
- Sleep Quality
- Fitness
- Smart Vibration
- Bluetooth
- APP
- Cloud

# O2 Vibe™ Sleep & Fitness Monitor

## Innovative Wearable Design

- The lightest and smallest SpO2 sleep monitor on the market  
Free Your Finger with the patented comfortable thumb sensor. No other pulse oximeter provides this capability!
- Dual Mode Sleep and Fitness Monitoring for day and night tracking and trending
- Better understand sleep quality using O2 Vibe Score
- Track how changes in fitness and lifestyle impact sleep quality
- Improve quality of sleep using customizable O2 thresholds and smart alarm settings
- Stay on track with daily step goals and smart alert reminders
- Track SpO2 while lightly exercising to maximize fitness regimens



## Smart Alert Vibrations Provide Sleep and Fitness Alerts

- Adjustable vibration setting gently wakes you when SpO2 falls below set threshold
- Don't interrupt partner during sleep with quiet personal smart vibrations
- Remind yourself to keep moving throughout the day to meet your step goals



## Mobile App and Secure Cloud Service

- Data and Trends can be synched with Secure Cloud Service automatically. Data and trends can be transferred automatically to trusted advisor or family member.
- Secure Cloud Service complies with all HIPAA regulations for privacy and protection of your data.

### Specifications:

Weight	35g (main unit)
Size	44x25x15mm (main unit)
Display	OLED
Battery	Rechargeable Lithium-polymer
Charge Time	Less than 2 hours to 90%
Wireless	Bluetooth 4.0 BLE
SpO2 Range	70% to 100%
Heart Rate Range	30 to 250 bpm
Vibrator	Triggered by low oxygenation or steps
Recorded Parameters	SpO2, heart rate, steps, motion
Mobile App for IOS	IOS 9.0 or above, iPhone 4s/iPad 3 or above
Mobile App for Android	Android 4.4 or above with Bluetooth 4.0 BLE

## Android Device Compatibility List

- Amazon Kindle Fire Tablet
- HTC One (M7)
- HTC One (M8)
- Nexus 5 and later
- Nexus 7 Tablet (2013 Version)
- Nexus 9 Tablet
- Motorola Moto X
- Samsung Galaxy Note 3 or later
- Samsung S4 or later
- Sony Ericsson Xperia Z2
- Sony Ericsson Xperia ZL
- Samsung S5 or later
- Samsung Galaxy Note 4 or later
- Motorola Moto X
- Huawei Honor 8, 5X
- Huawei P6 or later
- Huawei Mate8 or later

It is recommended to use O2 VIBE App with a validated device to ensure maximum reliability.

**O2 Vibe Monitor available online at for \$199 with free patient portal and dashboard.**