Chart #1

STETHOSCOPE USED WITH DIGITAL TRANSMISSION SOFTWARE

TELEMEDICINE STETHOSCOPE FEATURES COMPARISON	RNK PCP-USB w/PCP-SSP Software	3M™ Littmann® 3200 w/TeleSteth™ System (End-of-Life 1-2020)	Eko Core w/Al Screening App & Livestream Software	Eko Duo ECG & Digital Stethoscope w/AI Screening App & Livestream Software	JedMed OMNISteth™ w/MedImaging Software
In-Band* or Out-of-Band**	In-Band	Out-of-Band	Out-of-Band	Out-of-Band	Out-of-Band
Connection to PC	USB Port	Bluetooth	Bluetooth 4.0 LE	Bluetooth 4.0 LE	3.5mm mic port & BT 4
Power Source	USB Port Always On	AA Batteries (50-60 hours)	Rechargeable Li-Ion Battery (9 hours)	Wireless Induction Charging & Rechargeable Li-Ion Battery (9 hrs)	AAA batteries (8 hours)
Frequency Range & Filter Options	20Hz-2000Hz 5 Filter Settings (Software)	20Hz-2000Hz 3 Filter Settings (Manual)	20Hz-2000Hz (Bell/Diaphragm Double Side Head)	20Hz-2000Hz 4 Filter Settings - B/D/MR/EXT	20Hz-2000Hz 3 Filter Settings B/D/W (Manual)
Amplification***	Yes - To Max of PC's Mic & Speaker Amplification	24X Acoustic***	40X Acoustic***	60X Acoustic***	24X Acoustic***
Volume	PC Volume Control	External Manual	External Manual	External Manual Controls - 12	External Manual
Adjustment	Continuous Settings	Controls - 8 Settings	Controls - 7 Settings	Settings	Controls # Settings N/A
Pediatric Use	Yes - 1 1/3" (33mm) Raised Center for Pedi & Infant Use	Yes - Adult, Pedi and Infant Use	Yes - Pediatric Diaphragm Included – Switchable with Bell	Yes - Adult and Pedi - Not For Infant Use Under 10kg	Yes – 25 mm Pedi Head /40 mm Adult Head
Record & Forward	Yes – Separate Software	Yes – 12 exams	Yes – Record w/ App	Yes – Record w/ App	Yes - 160 exams****
Visualization	Yes w/Third Party Apps	Yes w/StethAssist™	Yes – On App	Yes – On App & Web Dashboard	1.5 in screen
OS Compatibility	Windows	Windows	Android/IOS/Win	Android/IOS/Win	Info Not Available
API/SDK	Win/WebRTC- Chrome/Android/IOS	Win XP-7/Fedora Core/Android	Win/IOS/Android	Win/IOS/Android	Info Not Available
AI Decision Support	No	No	Murmur/AFib Screening	Murmur/AFib Screening	No

*In-Band Software transmits the audio signal through the audio channel of the video conference platform

**Out-of-Band Software transmits the audio signal on a separate channel from the video conference platform

***Amplification ratings for stethoscopes are shown as maximum amplification vs traditional acoustic stethoscopes. All digital stethoscopes have adjustable volume levels and maximum volumes will vary.

****Organ button selects body area templates with preset frequency ranges for recording.

Disclaimer: Chart data is based on currently available information and may be subject to change with product updates or pricing changes

www.mobildrtech.com

www.telemedicinesupply.com

Chart # 2

COST COMPARISON - STETHOSCOPES USED WITH DIGITAL TRANSMISSION SOFTWARE

TELEMEDICINE STETHOSCOPES – TOTAL COST OF OWNERSHIP	RNK PCP-USB w/ PCP-SSP Software	3M [™] Littmann [®] 3200 w/ TeleSteth [™] System (Software End- of-Life 1-2020)	Eko Core w/ Al Screening App & Livestream Software Service	Eko Duo ECG & Digital Stethoscope w/ Al Screening App & Livestream Software Service	JedMed OMNISteth™ w/ MedImaging Software
Stethoscope	\$489	\$399 X 2 = \$798*	\$249	\$349	\$480 X 2 = \$960**
Software License	\$96 one-time / stethoscope	\$750 X 2 = \$1500/yr	None	None	\$1100 lifetime license
Recurring Fees	None	None	Store and Forward Free w/Livestream & Murmur/AFib \$1200/yr w/Video, Livestream & Murmur/AFib AI Screen \$2400/yr	Store and Forward Free w/Livestream & Murmur/AFib \$1200/yr w/Video, Livestream & Murmur/AFib AI Screen \$2400/yr	\$50/yr/stethoscope
1 Year Cost	\$585	\$2298	\$249 - \$2400	\$349 - \$2400	\$2160
3 Year Cost	\$585	\$5298	\$249 - \$4800	\$349 - \$4800	\$2360
5 Year Cost	\$585	\$8298	\$249 - \$120 0 0	\$349 - \$12000	\$2560

* Stethoscopes required at sending and receiving sites for real-time auscultation

** Second stethoscope for JedMed not required but recommended by manufacturer for best quality sound

 Disclaimer: Chart data is based on currently available information and may be subject to change with product updates or pricing changes

 www.mobildrtech.com
 www.telemedicinesupply.com

 support@mobildrtech.com

Chart # 3

STETHOSCOPE USED AS EXTERNAL MICROPHONE W/O SOFTWARE

TELEMEDICINE STETHOSCOPE FEATURES COMPARISON	RNK PCP-USB w/o Software	RNK PCP-1 w/o Software	Cardionics E-Scope II w/o Software	Thinklabs® Digital One w/o Software
Microphone Switching	Switch Mic in VC Software	Switch Mic in VC Software	Switch Mic in VC Software	Switch Mic in VC Software
Connection to PC	USB Port	3.5 mm Mic Port	1/8 in male Stereo Jack & USB Accessory Connector	3.5 mm Mic Port
Power Source	USB Port Always On	Mic Port Always On	AAA Battery (4.5 mo Regular Use)–Auto Off After 1.5-2 Min	Rechargeable Li-Ion Battery (100-125 exams)-Auto Off Can Be Set To 1-10 Minutes
Frequency Range & Filter	20Hz-2000Hz	20Hz-2000Hz	20Hz-2000Hz	20Hz-2000Hz
Options	No Frequency Selections w/o Software	No Frequency Selections w/o Software	Heart and Breath Filter Settings (Manual)	5 Filter Settings (Manual)
Amplification*	Yes – To Max of PC's Mic & Speaker Amplification	Yes – To Max of PC's Mic & Speaker Amplification	30X Acoustic*	100X Acoustic*
Volume Adjustment	PC Volume Control Continuous Settings	PC Volume Control Continuous Settings	External Manual Controls 64 Settings	External Manual Controls 10 Settings
Pediatric Use	Yes – 1 1/3" (33 mm) Raised Center For Pediatric & Infant Use	Yes – 1 1/3" Raised Center For Pediatric & Infant Use	Supplied w/ Specialist Adult Diaphragm That May Be Changed to Specialist Adult or Specialist Pediatric Bell	46 mm Adult Diaphragm Has Been Used for Pediatric Patients
Record & Forward	Yes – Separate Software Required	Yes – Separate Software Required	Yes w/Connecting Cable (cat. No. 711-7128)	Yes w/Thinklink Connector
Visualization	No –	No –	No –	No –
Ease of Use**	Yes w/ Third Party Apps	Yes w/ Third Party Apps	Yes w/ Third Party Apps 2	Yes w/ Third Party Apps
Cost	\$489	\$399	\$335	3 \$499
CUSI	ر ن ب ر	3333	,,,,,	ς τ τς

*Amplification ratings for stethoscopes are shown as maximum amplification vs traditional acoustic stethoscopes.

All digital stethoscopes have adjustable volume levels and maximum volumes will vary.

**Ease of Use Scale 1=Easiest-No External Buttons or Switches 2=External Manual Controls on Stethoscope Head 3=Most Technical for End User