

# USB ExamCam 2.0

The USB ExamCam 2.0 turns a high resolution, auto-focus web camera into a high quality, general purpose exam camera at an unprecedented price point. Using custom designed parts from RNK Products, a leading manufacturer of stethoscopes and medical devices, MobilDrTech has assembled a functional, high quality general examination camera.

The USB ExamCam 2.0 is available with free software that will provide picture snap and video recording capabilities. It has a circular LED light ring attached to the camera for high quality viewing.



# Features

- **Light Ring** – Permanently affixed circular light ring with 6 LEDs for clear, crisp viewing of skin or throat



See mp4 video clip at: <https://cdn.shopify.com/s/files/1/1693/8993/files/Throat.mp4?74>

- **Easy Grip Detachable Handle** – Attaches to webcam for one handed operation and always correct camera orientation – Fitted with ring for hanging camera on a hook when not in use – Ring is removable to expose a 1/4"-20 thread hole at the bottom that allows mounting to any tripod or monopod for fixed position use.
- **Single Cable** – USB cables from the base camera and the light ring are constructed into a single cable and connect through a 2-port USB Hub device to the PC or Laptop with a single USB plug. The USB ExamCam 2.0 ships with 9 ft USB extension cable for an effective cable length of 15 feet. This insures ease of access to the patient in office and clinic settings.
- **Say Ahh Band Attachment** – Designed to fit snugly around the light ring so that a tongue depressor can be inserted under the camera – This makes one handed operation possible when conducting a throat exam.
- **Free Software Available** – Adjustable focal length with picture snap and video clip capture.
- **USB Plug-and-Play** – The camera may be used standalone or with software.
- **LED Lighting** – Powered by USB port on computer – Comes on automatically when the camera is plugged into the USB port – Lighting may be unplugged at 2-port USB hub if only ambient room light is desired.



MobilDrTech is an approved Value Added Reseller / Distributor for products from RNK. All products shown are available on the MobilDrTech website at [www.telemedicinesupply.com](http://www.telemedicinesupply.com) or contact us for more information at the email/phone shown below.

## Pricing

**The USB ExamCam 2.0 is a fully assembled and finished general examination camera. Partial component parts kits are available for customers interested in DIY. Just let us know what works best for your program.**

**Fully Finished Exam Cam – Plug-and-Play – Ready to Use for \$249**

## Contact Us

MobilDrTech, Inc.  
14090 SW Fwy, Ste 300  
Sugar Land, TX 77478  
281-340-2013 Office  
281-340-2001 Fax  
[sales@mobildrtech.com](mailto:sales@mobildrtech.com)  
[www.mobildrtech.com](http://www.mobildrtech.com)

## Products

[www.usbsteth.com](http://www.usbsteth.com)  
[www.usbexamcam.com](http://www.usbexamcam.com)  
[www.telemedicinesupply.com](http://www.telemedicinesupply.com)