

VitalsRx Monitor


Settings for English and Spanish languages

Finally, a pocket-size, professional quality Vitals Monitor that is reasonably priced, FDA approved and does not require the use of specific providers or subscription service. Data can be shared with any physician, caregiver or family member from a free App on IOS or Android devices.



Remote Physiologic Monitoring (RPM) and Chronic Care Management (CCM) are two rapidly growing areas of medicine. In early 2018, CMS "unbundled" the CPT code for Remote Patient Care, 99091. This

separated qualifying services from restrictive telehealth and telemedicine reimbursement guidelines. It has opened the door for healthcare providers to be reimbursed for some RPM services including data collection and care coordination for patients, including those at home.



HEALTH FEATURES & METRICS


BodiMetrics™ VitalsRx captures vital health metrics in the palm of the hand in less than 30 seconds! Vitals (ECG - Heart Rate - Respiration - Heart Rate Variability - Blood Oxygenation - Temperature); pedometer and medication reminders combined with algorithms and indices provide the information needed to support a healthy lifestyle, chronic care management or remote patient monitoring programs.

- 20 Second ECG Waveform
- Continuous ECG to 24 Hrs
- 20 Second SpO2 Capture
- Continuous SpO2 to 15 Hrs
- Overnight Sleep Screening
- Heart Rate & Variability
- 3 Second Thermometer
- Stress Relief Bio-Feedback
- Pedometer Steps
- Medication Reminders

USE CASES – RPM/CCM

- Sleep Apnea Screening
- Cardiac Screening
- Cardiac Recovery
- Post Discharge Care
- Senior Care & Wellness
- 24 Hour Holter Monitoring
- Pharmaceutical Research
- Medication Adherence
- Behavioral Health
- Employee Health

BodiMetrics™ VitalsRx tracks and reports on continuous ECG (24 hour Holter) and Blood Oxygenation (SpO2 for up to 15 hours) providing effective monitoring and screening for heart irregularities and sleep related health issues like sleep apnea. Optional Bluetooth Blood Pressure wrist cuff provides the user seamless Blood Pressure metrics and syncing to VitalsRx mobile app and secure cloud.



Telemedicine...
in the palm of your hand.



BodiMetrics™ Cloud is a powerful system that provides patients and clinicians secure and easy connections. Data can flow from one person to another seamlessly and freely under authorization. Both patients and clinicians benefit from the availability of data.

- Data can be synced to cloud automatically
- Data can be shared with doctor automatically
- Data can be stored in cloud over long term for health tracking
- Data is encrypted to protect user privacy




VitalsRx (4+)
BodiMetrics, LLC
Free



Patient Monitoring Dashboard – Free Cloud Service Available with All Monitors – Follow Unlimited Number of Patients from Multiple Caregivers:

Performance Monitor™ and VitalsRx™ Health Monitors

Patient Portal



Monitor Tab – Admin & Viewer

Bodimetrics										Admin: mgoetting@mindstrel.com
Monitor										Logout
Dashboard										
PWI View										
Accounts										
Download										
Settings										
Welcome!										
Most Recent Observations										
Show <input type="text" value="20"/> entries										Search: <input type="text"/>
Last Name	First Name	Medical ID	Temp	SBP/DBP	Pulse	Resp	SpO2	PWI	Alerts	
AL	ANDROID	987654	97.3	--	77	10	99	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
SMITH	BILL	12345	97.2	109	66	8	97	59	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
SMITH	BILL	12345	97.5	-	-	-	-	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
ADAMS	JACK	1009999	97.2	98	62	10	98	10	<input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
WILSON	JACK	1001111	98.1	117	68	13	99	90	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
	JACK	1009999	97.7	120	83	9	99	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
FDA	JANE	9876543210	-	-	69	-	96	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
EYRE	JANE	88888888	97.5	132	82	10	99	19	<input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
	JANE	88888888	97.7	146	81	12	97	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
	JOCELYN	1	97.9	112	77	11	99	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
WILSON	JOE	112233445566778899	97.5	140	100	11	96	28	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
yahoo	Karen	123	-	-	-	-	-	-	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	

VitalsRx Monitor available online for \$499 with free patient portal and dashboard.