

Superior ECG Technology for the 21st Century

Anyone, Anytime, Anywhere



Finally, a state-of-the-art 12-lead ECG with 21st century technology! With continuing software and diagnostics updates and unlimited cloud storage, you will always be up-to-date for years to come!

PCA 500 Platform – The world's first 12-lead ECG that enables anyone to complete a hospital-quality ECG, even without any training.


- ▶ **PCA 500** – Compact & portable.
- ▶ **Patented Prepositioned Electrode Strip** – Hassle-free placement – anyone can do it!
- ▶ **QT ECG App** – Record a high quality 12-lead ECG at the touch of a button, on a smartphone or tablet.
- ▶ **QT Cloud** – Unlimited HIPAA & GDPR compliant cloud storage and database management.
- ▶ **QT Diagnostics** – Instant preliminary computer interpretation using the most up-to-date algorithm.




- Weight 67 gm (2 oz)
- 1,000 Hz sampling rate
- Rechargeable Li-ion battery
- Splash-proof



- Instant preliminary computer interpretation
- Customizable notifications to physicians
- Easy-to-use cloud dashboard
- Access all ECG records remotely



- Simple electrode placement by anyone with minimal training
- AHA and IEC compliance for clinical standard electrode placement
- Four electrode strip sizes for adults to ensure correct fit
- Sizing guide to select the right size electrode strip



- Mobile apps for iOS and Android phones and tablets
- Intuitive user interface with step-by-step guidance
- Lead-off detection and alert
- ECG report in pdf file for printing or sending by email