

Biological Indicators & Real Time Automatic Reader

**ULTRA
RAPID
20 MIN.**

Biological Indicators (Ultra & Super rapid)
and Incubators

Designed for Hospitals, ASC's, Physicians & Dental Practices

With increasing demands for faster sterile instruments
turnaround in healthcare facilities, sterilization results need
to be available as soon as possible to reliably verify the sterility
of the load.



Tuttnauer
Innovation · Legacy · Partnership

Biological Indicators & Real Time Automatic Reader

Automatic reading and traceability system for sterilization

- 20 min or 1 hour process times
- Multiple process times can be used simultaneously
- Ease of use
- Thermal printer
- USB Connection with free software download
- One BI & reader can be used in both gravity and pre/post vacuum sterilizers



Biological Indicators

WTL198-0058



1 HOUR

GRAVITY/VACUUM
RAPID

WTL198-0075



1 HOUR

GRAVITY/VACUUM

WTL198-0072



20 MIN

GRAVITY/VACUUM
ULTRA RAPID

WTL198-0077



20 MIN

GRAVITY/VACUUM


Chemical Indicators and Tape

WTL198-0081



MULTI-VARIABLE INDICATOR TYPE 4

WTL198-0082



INTEGRATORS TYPE 5

WTL198-0053



SELF-ADHESIVE TAPE TYPE 1

Tuttnauer
Innovation · Legacy · Partnership

Learn from our experts.

Visit our blog: tuttnauer.com/blog



Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720

E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

