

New Release

2540MKA –

The First Manual Autoclave with Closed Door Drying



Our most recent development, the 2540MKA, is the first manual autoclave in the market featuring **closed door active drying** at an affordable price. The 2540MKA sterilizes and dries your instruments in a single action, saving you the trouble of supervising the cycle.

The 2540MKA features heater auto-shutoff at completion of sterilization cycle, saving energy and increasing your safety and productivity. Before you leave the clinic, set the 2540MKA to Night Mode, start a cycle and the instruments will be sterile and ready to use when you return.

Benefits:

- Closed door active drying air pump
- Single Action Sterilization Cycle
- Automatic night mode shutoff time and energy saver

Safety Features:

- Steam pressure activated lock prevents door from opening if chamber is pressurized
- Double protection mechanism prevents door from opening at high pressure
- Automatic safety shutoff prevents overheating
- Insufficient water detector with automatic shutoff
- Overheat Safety protection
- 0.2µm HEPA air filter provides sterile, bacteria-free air for drying

Complies with strictest international standards and directives.



2540MKA Panel Features:



- ① Cycle Status Indicator
- ② Pressure Indicator
- ③ Power Switch
- ④ Sterilization Timer
- ⑤ Water Pump Switch
- ⑥ Drying and Sterilization Timer
- ⑦ Night Mode Switch
- ⑧ Temperature Selector

2540MKA Technical Specifications:

Chamber Dimensions \varnothing x L (mm)	250x467
Chamber Volume	23L / 6.1gal
Cycle Time	20 Min.
Standard Cassettes Capacity	3 Full
No. of Trays	4
Tray Dimensions WxHxD (mm)	168x20x414
Voltage(V) / Freq.(Hz)	230V / (50/60 Hz)
Current(A) / Power(W)	10.0A / 2200W
Overall Dimensions WxHxD (mm)	508x362x550
Autoclave Weight (Kg.)	28 kg



Tuttnauer Europe b.v.

Paardeweide 36,4824 EH, P.O.B. 7191, 4800 GD Breda, The Netherlands
 Tel:+31(0)76 5423510, Fax:+31(0)76 5423540, E-mail: info@tuttnauer.com, www.tuttnauer.com