

Sensitivität / Sensitivity

AT-Nr. / AT-No.	AT-Nr. / AT-No.	Selbsttest / self- test	Ref.Nr./ ID-No. *	Hersteller / Manufacturer	Testname / Test name	Zielantigen / target antigen	Cq ≤25	Cq 25-30	Cq ≥30	Gesamt- Sensitivität / total sensitivity
AT883/21			CO-05	Jiangsu Weil Biotech Co., Ltd.	COVID-19 Schnelltestgerät	N	100,0%	95,0%	40,0%	86,0%
AT770/21			n.a.	Jiangxi Province Jinhuan Medical Instrument Co., LTD.	COVID-19 antigen Schnelltest (Kolloidales Gold) / Novel Coronavirus(SARS-CoV-2) Antigen Rapid Detection Kit	N	100,0%	95,0%	60,0%	90,0%
AT882/21			n.a.	Jiaxing Wisetest Bio-Tech Co., Ltd.	Saliva SARS-CoV-2 (2019-nCoV) Antigen Test Kit (Nanocarbon Assay)	N	85,0%	15,0%	10,0%	42,0%
AT559/21			n.a.	Jinan Babio Biotechnology Co., Ltd.	SARS-CoV-2 Antigen Rapid Detection Kit (Colloidal Gold Method)	N	100,0%	35,0%	0,0%	54,0%
AT286/20	AT127/21	n.a.	n.a.	Joinstar Biomedical Technology Co., Ltd (CIV care impuls Vertrieb)	COVID-19 Antigen Schnelltest (Colloidal Gold)	N	100,0%	60,9%	0,0%	64,0%
AT527/20			COV-AG-20 / G10313	JOYSBIO (Tianjin) Biotechnology Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) SPUCKTEST	N	100,0%	90,0%	0,0%	76,0%
AT692/21			G10322	JOYSBIO (Tianjin) Biotechnology Co., Ltd.	Joyshio SARS-CoV-2 Antigen rapid test kit -PRO (Colloidal Gold)	N	100,0%	35,0%	0,0%	54,0%
AT127/20	AT1284/21 / AT1159/21		B63000	Labnovation Technologies, Inc.	Labnovation SARS-CoV-2 Antigen Rapid Test Kit (Immunochromatography)	N	94,1%	17,4%	0,0%	40,0%
AT1054/21			n.a.	Lifecorn Biotech Limited	COVID-19 Rapid Test Cassette Antigen Test Kit	N	100,0%	100,0%	100,0%	100,0%
AT145/20			PR210101	Lumigenex (Suzhou) Co., Ltd.	PocPoc SARS-CoV-2, Antigen Schnelltest Set (Kolloidales Gold)	N	100,0%	65,2%	0,0%	64,0%
AT151/20			71110	LumiQuick Diagnostics, Inc.	Quick Profile Covid-19 Antigen Test Card	N	100,0%	91,3%	20,0%	80,0%
AT075/20			SPEC-32312 R2 ART-00571 R4	LumiradX	LumiradX SARS-CoV-2 Ag Test	N	100,0%	52,2%	0,0%	60,0%
AT012/20			76748	MEDsan GmbH	MEDsan® SARS-CoV-2 Antigen Rapid Test	N	100,0%	47,8%	0,0%	58,0%
AT214/20			CP01810011	Merlin Biomedical (Xiamen) Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Cassette	N	100,0%	82,6%	0,0%	72,0%
AT790/21			CP02150011	Merlin Biomedical (Xiamen) Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Cassette	N	100,0%	55,0%	0,0%	62,0%
AT047/20			0230005SP	Molab GmbH	mfo-screen Corona Antigen Test	N	100,0%	47,8%	0,0%	58,0%
AT198/20			07AG6020B	MP Biomedicals Germany GmbH	Rapid SARS-CoV-2 Antigen Test Card	N	100,0%	43,5%	0,0%	54,0%
AT823/21	AT1281/21		0000000168	MSP bodmann GmbH	MSP SARS-CoV-2 Antigen-Schnelltest	N	100,0%	20,0%	0,0%	48,0%
AT021/20	AT1293/21		243103N-20	nal von minden gmbh	NADAL® COVID-19 Ag Schnelltest	N	83,3%	13,0%	0,0%	36,0%
AT593/21			500200	Nanjing Liming Bio-Products Co., Ltd.	StrongStep® SARS-COV2 Antigen Rapid Test	N	95,0%	80,0%	0,0%	70,0%
AT462/20			n.a.	Nanjing Norman Biological Technology Co., Ltd	Novel Coronavirus (2019-nCoV) Antigen Testing Kit (Colloidal Gold)	N	94,1%	26,1%	0,0%	44,0%
AT678/21			RQ005	Nanjing Syntogene Medical Technology Co., Ltd.	Syntogene Neues Coronavirus (SARS-CoV-2) Antigen-Schnelltestmethode (Collate Gold-Methode)	N	100,0%	35,0%	0,0%	54,0%