



- Lightweight, easy-to-handle and set up
- Folds in half, carries like a suitcase*
- Durable welded construction
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included



| Item Number | Length | Width | Folded Size | Weight | Capacity** |
|-------------|--------|-------|-------------------------------|---------|------------|
| SFW230* | 2' | 30" | 24" x 16" x 4 $\frac{7}{8}$ " | 10 lbs. | 800 lbs. |
| SFW330 | 3' | 30" | 36" x 16" x 4 $\frac{7}{8}$ " | 15 lbs. | 800 lbs. |
| SFW430 | 4' | 30" | 48" x 16" x 4 $\frac{7}{8}$ " | 20 lbs. | 800 lbs. |
| SFW530 | 5' | 30" | 60" x 16" x 4 $\frac{7}{8}$ " | 24 lbs. | 800 lbs. |

Not sure which ramp length you need? Refer to pages 18 - 19 to help you.

*Handle not included on SFW230. ** Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles).

LIMITED LIFETIME WARRANTY

MOBILITYCARE
Since 1992

PVI

Phone: (03) 9568 8383

39 Cambro Road, Clayton VIC 3168

info@mobilitycare.net.au

www.mobilitycare.net.au

(800) 554.7267 visit us online

com