



EICKEMEYER®  
I.V. Fluid Warmer

veterinary technology for life  
**Eickemeyer**®

# Nice and warm



## **EICKEMEYER® IV Fluid Warmer**

The EICKEMEYER® blood and infusion warmer reduces the temperature difference between IV fluid and the patient and minimises the occurrence of hypothermia.

- ▶ For warming up all intravenously administered liquids
- ▶ Easy to handle
- ▶ LED display with temperature display
- ▶ Precise temperature control between 35 °C and 42 °C (0.1 °C steps)
- ▶ Alarms: overheating, temperature too low, system error, door open, heating reminder
- ▶ Two fastening options: with a lanyard or fastening clamp
- ▶ Warm-up time: less than 2 min
- ▶ Power supply: 100 – 240 V
- ▶ Mains voltage: 50/60 Hz
- ▶ Power: 85 VA
- ▶ Dimensions (in cm): L 18 x W 7 x D 3.5
- ▶ Weight: 200 g (without fastening clamp), 400 g (with fastening clamp)

**402092**

### **EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**