

DKT.P5K01.5000

Filter: P5K01

Veterinarian Laser Safety Doggles® Eyewear Kit

Color: Dark Green
Material: Polycarbonate
VLT: 10%




Product is designed to protect animal's eyes from laser exposure during laser therapy. Kit includes three sizes to fit most pets from a small cat to a large dog.

Optical Density (OD):
 OD 5+ @ 800 - 1100 nm

- Kit includes D01, D02, and D03 in soft pouches.
- Eyewear is sold individually or as a kit.
- Filter cover most therapeutic lasers.

Doggles® is a registered trademark of Doggles®, LLC.




Visual interpretation of data. Please refer to the OD specifications for the actual protection level.

Frame Specifications

Sold individual or as a kit

Part No. & Dimensions	D01. P5K01.5000	D02. P5K01.5000	D03. P5K01.5000
Pet Size	xSmall Pets Below 10 lbs.	Small/Medium Pets 9-25 lbs.	Large Pets 50-95 lbs.
a. Lens Width	31 mm	47 mm	60 mm
b. Lens Height	22 mm	31 mm	38 mm
c. Between Lenses	16 mm	23 mm	26 mm
d. Temple Length	Adjustable	Adjustable	Adjustable
e. Nylon Strap	Adjustable	Adjustable	Adjustable
f. Frame Height	25 mm	44 mm	53 mm
g. Frame Length	Adjustable	Adjustable	Adjustable
h. Frame Width	93 mm	130 mm	160 mm



DKT.P5K01 Box



Box size - 7" x 6" x 2.25"



D03.P5K01.5000



Authorized Laservision Signature

Date