

EXTRACTION FORCEPS




Extraction forceps are used to further loosen teeth with rotational and tilting movements and finally to extract the teeth in combination with a fulcrum.

The **Extraction Forceps No. 1** is particularly helpful for mandibular cheek teeth extractions but can be used for the upper jaw as well.

The **Extraction Forceps No. 2** is specifically designed for the wider maxillary cheek teeth.




The Extraction Forceps can be disassembled and reassembled again in many different ways like the spreaders.

For horses with small heads we developed the "**Pony Model**". The actual version is dedicated to **Professor Padraic Dixon** from the University of Edinburgh who, at the end of the 20th century, made great efforts to re-introduce oral extraction techniques.



REVERSE FULCRUM

Thanks to this NEW Reverse Fulcrum Extension there is no extra Reverse Fulcrum Forceps needed. This system works both for lower and upper cheek teeth. It is especially recommended for 06 and 07's. It can be a quite helpful approach for lower 09 and 10 as well due to the curved nature of these teeth.



Art.-No. #	Description
PEGA_321020	Molar Forceps – Pony
PEGA_321001	Molar Forceps No. 1 – large
PEGA_321002	Molar Forceps No. 2 – large
PEGA_321010	Reverse Fulcrum Extension (fits Molar Forceps No. 1 and No. 2)