



TECHNICAL DATA SHEET

REPROJELLY

Product description

ReproJelly is a sperm-safe lubricating gel.

Application

It is recommended for AI, rectal palpation and lubrication of AVs in animal species.

Properties

- Composition: purified water, hydroxyethyl cellulose, bacteriostatic and fungicidal agents, stabilizers
- Appearance: clear to light yellow gel
- Odor: neutral
- Produced under GMP norms; components are certified corresponding to Ph EUR
- No ingredients of animal origin
- Iso-osmotic
- Tested for sperm compatibility
- Viscosity optimized to coat gloves, AI equipment and instruments

Storage

+5 to +25°C

Order information

ReproJelly, 100 ml tube	REF. :	11907/0100
ReproJelly, 3 l canister	REF. :	11907/3000
Dispenser for 3 l canister	REF. :	11907/3100

