
## REVERSE FULCRUM

Thanks to this NEW Reverse Fulcrum Extension there is no extra Reverse Fulcrum Forceps needed. This system works both for lower and upper cheek teeth. It is especially recommended for 06 and 07's. It can be a quite helpful approach for lower 09 and 10 as well due to the curved nature of these teeth.





ArtNo. #	Description
PEGA_321020	Molar Forceps – Pony
PEGA_321001	Molar Forceps No. 1 – large
PEGA_321002	Molar Forceps No. 2 – large
PEGA_321010	Reverse Fulcrum Extension (fits Molar Forceps No. 1 and No. 2)