SYS-70 Infusion Pump The champion of single pump

Innovative and Safe



Ingress protection by IP34

Adapt to the needs of clinical use in various environments



Fits different applications



Applies to different medications and treatment methods



Large capacity Li-ion battery

Supports up to 20 hours usage



Dual occlusion alarm system

Upstream and downstream occlusion alarm with 13 levels



Intermittent air bubbles alarm

Detects air bubbles with great sensitivity at 7 levels



Information Interconnection

SYS-50/52and SYS-70 respectively have an optional built-in WiFi module. It may connect via WLAN or LAN connection to any Hospital Information System (HIS) a They are connected to Medcaptain's MP-900 infusion central monitoring system to implement centralized monitoring of all infusions. All SYS series pumps can also access the Clinical Information





MEDCAPTAIN[]

MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.



Add: 12th Floor, Baiwang Research Building, No.5158 Shahe Villi, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC Tel: +86 755 26953369 E-mail: info@medcaptain.com
Fax: +86 755 26001651 Website: www.medcaptain.com

MEDCAPTAIN is a trade mark of Medcaptain Medical Technology Co., Ltd.

Specifications are subject to change without prior notice. Add: 12th Floor, Baiwang Research Building, No.5158 Shahe West Road, Xili, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA



Distributor:

SYS-50 / SYS-52 Syringe Pump **SYS-70 Infusion Pump**

MEDCAPTAIN[]

SYS-50/52 Syringe Pump

Simple



Touch what you see. See the parameters you need on the screen



A standard fitted side clamp will prevent infusion tube entanglement and enhance infusion overview and control



Continuous titration

Changing the flowrate during infusion without stopping the infusion



Open design

Facilitating convenient and fast syringe installation and constant observation of remaining liquid in the syringe





Ingress protection by IP34

Meets different kinds of environmental requirements



Drive Head Protection Design

Which is easy for transportation and can avoid bolus if pumps fall down



Large capacity Li-ion battery

Supports up to 20 hours usage





-దౖ- Pre-Occlusion Alarm System

Which can warn clinical personnel before occlusion alarms actually happen to ensure safe infusion

Closer to technology, closer to you **Smart**

- Different infusion modes integrated with SYS-50/52 and SYS-70 fulfil clinical demands
- Connection to Center Monitoring System enables real time monitoring and centralized infusion management
- Docking Station makes one charging cable possible for all pumps





SYS-52 Syringe Pump





EN 1789/standard for ambulance

SYS-50/52 syringe pump, SYS-70 infusion pump designed for ambulance and medical transportation, during repeated tests, they acquired EN 1789 certification, which guarantees stable infusion during transportation.

SYS-70 Infusion Pump The champion of single pump

Simple and Significant



3.5" Touch screen operation

Touch what you see. See the parameters you need on the screen



Automatic pump door & Automatic anti free-flow system

One touch to open the door and clamp. Convenient and easy operation



Flexible and lightweight



- 01 Upstream Occlusion Alarm
- 02 Innovative Pressure-Adjustable module
- 03 Intermittent Air Bubbles Detector
- 04 Downstream Occlusion Alarm