

Specialized V-Top PRO ErgO Laparoscopic surgical table



- Specialized V-Top PRO ErgO laparoscopic table
Body: pulvercoated metal frame
Table top : V-Top PRO ErgO operationtable 133x52cm
- * Electrical trendelburg and antitrendelburg tilt
 - * Electrical lateral tilt
 - * Electrical height adjustment :
 - 65-105cm without wheels
 - 75-115cm with wheels


4100 Specialized V-Top PRO ErgO operatietafel

Optional accessories: nr 1,3,4,6,7,8,9,11,12.


V-Top PRO ErgO accessories




01 Dentalextension



02 Draincontainer



03 extension parts




04 Tubing holder




05 Fixation frame



06 RVS bucket



07 Fixation clips



08 Wheel set




09 Rubbermat



10 Built-in weighing scale



11 Cable protector



12 Infusion stand