



Versa LED High Speed Handpiece with Swivel



Item #54159

KEY FEATURES:

CERAMIC BALL BEARINGS

- Superior torque and long lifespan
- Specially designed for smooth and quiet operation

TRIPLE WATER SPRAY

Effectively cools the entire surface of the bur to prevent overheating, resulting in enhanced patient comfort.

DIRECT LED

The LED is optimally positioned on the face of the handpiece to generate a larger light pattern than a traditional fiber-optic system.

TITANIUM SURFACE

Solid titanium construction that is built to last. Microtextured and scratch-proof surface reduces the possibility of slippage.

GENERATOR

An innovative micro-generator producing electricity for the LED is built into the handpiece. This eliminates the need for a fiber-optic system.

COMPATIBILITY

Compatible with a wide variety of coupling systems and adapts well to most existing treatment units.

SPECIFICATIONS:

HEAD SIZE (W x H) MM	10.8 x 13.0
POWER (WATT)	23
DRIVE AIR PRESSURE (PSI)	40-60
SPEED (RPM)	410,000
WATER SPRAY	Triple Ports
COLOR TEMPERATURE (K)	5,500
BRIGHTNESS (LUX)	25,500
O / MAX. LENGTH (MM)	1.6 / 25
CHUCK TYPE	Push Button
ROTARY INSTRUMENTS	FG burs according to EN ISO 1797-1

Item	Product Description	UOM
54159	Versa LED High Speed Handpiece with Swivel	Each

**NO LIGHT BULB,
NO FIBER-OPTIC TUBING,
JUST BRIGHT ILLUMINATION!**