













Woodcast Product Catalog 2021


Woodcast 1 mm Ribbon	REF	Finish	Size	Batch	Colors
For finger and thumb splinting and reinforcement. If you make circumferential casts, apply Ribbon stripes on top.	D4500	Non-vented	1mm x 4cm x 5m	2 pcs	
	D8750	Non-vented	1mm x 8cm x 7.5m	1 pcs	



Woodcast 2 mm	REF	Finish	Size	Batch	Colors
For all upper limb applications and smaller in lower limb applications.	C1434	Vented	2mm x 14.5cm x 34cm	20 pcs	
	C1480	Vented	2mm x 14.5cm x 80cm	10 pcs	
	C1434S	Non-vented	2mm x 14.5cm x 34cm	20 pcs	
	C1480S	Non-vented	2mm x 14.5cm x 80cm	10 pcs	



Woodcast 2 mm Soft	REF	Finish	Size	Batch	Colors
Ideal for hand therapy and focused rigidity casting together with Woodcast 2mm.	F1434	Vented	2mm x 14.5cm x 34cm	20 pcs	
	F1480	Vented	2mm x 14.5cm x 80cm	10 pcs	
	F1434S	Non-vented	2mm x 14.5cm x 34cm	20 pcs	
	F1480S	Non-vented	2mm x 14.5cm x 80cm	10 pcs	





Woodcast 2 mm Folded	REF	Finish	Size	Batch	Colors
Folded version of 80 cm long Woodcast 2mm that fits in Woodcast Express Heater	C1480F	Vented	2mm x 14.5cm x 80cm	10 pcs	



Woodcast 2 mm 4in1	REF	Finish	Size	Batch	Colors
Pre-cut splint for most common upper limb applications	C1441	Vented	2mm x 14,5cm x 29cm (pre-cut in shape)	20 pcs	











Make seven applications without cutting!


Woodcast OneStep	REF	Finish	Size	Batch	Colors
Pre-cut splint for ankle and elbow	C1465 (Elbow)	Vented	2mm x 14.5cm x 65cm	10 pcs	
	B1565 (Ankle)	Vented	4mm x 15cm x 65cm	5 pcs	




Heat and apply! No cutting needed

Woodcast 4 mm	REF	Finish	Size	Batch	Colors
Thick splint that provides a good support. Used for lower limb applications.	B1040	Vented	4mm x 10cm x 40cm	20 pcs	
	B1080	Vented	4mm x 10cm x 40cm	10 pcs	
	B1240	Vented	4mm x 12.5cm x 40cm	20 pcs	
	B1240S	Non-vented	4mm x 12.5cm x 40cm	20 pcs	 
	B1280	Vented	4mm x 12.5cm x 80cm	10 pcs	
	B1540	Vented	4mm x 15cm x 40cm	20 pcs	
	B1580	Vented	4mm x 15cm x 80cm	10 pcs	




Woodcast 4 mm Folded	REF	Finish	Size	Batch	Colors
Folded version of 80 cm long Woodcast 4mm that fits in Woodcast Express Heater	B1580F	Vented	4mm x 15cm x 80cm	8 pcs	




Stockinette for Woodcast	REF	Size	Batch	Colors
<p>Stockinette is used for protecting the skin. Use a minimum of 2 layers with Woodcast 2mm.</p> 	X1022	10cm x 22.8m	1 pc	White
	X1522	15cm x 22.8m	1 pc	White
	X2522	2.5cm x 22.8m	1 pc	White
	X5022	5cm x 22.8m	1 pc	White
	X7522	7.5cm x 22.8m	1 pc	White

Adhesive AirPadding	REF	Size	Batch	Colors
Breathable padding material. Use always with Woodcast 4mm.	X4350A	40cm x 3.5m	1 pc	Black



Extremely light and breathable

Edging tape for Woodcast	REF	Size	Batch	Colors
For protecting and padding of the edges.	X0315	3cm x 7.5m	2 pcs	Black



Hook for Woodcast	REF	Size	Batch	Colors
Adhesive Hook for orthoses that offers a unique, strong hold.	X5050	50mm x 22.8m	1 pc	Black



Stays in place with our unique polymer bond

Elastic loop for Woodcast	REF	Size	Batch	Colors
Elastic loop for making orthoses. Use together with Hook.	X0250	25mm x 7.5m	2 pcs	Black
	X5252	50mm x 16m	1 pc	



Loop for Woodcast	REF	Size	Batch	Colors
Loop for making orthoses. Use together with Hook.	X5251	50mm x 22.8m	1 pc	Black



Freeze spray	REF	Size	Batch	Colors
Used to cool down splint faster	X0400	400 ml	1 pc	-



Scissors	REF	Size	Batch	Note
Economy Scissor	MPX2020		1 pc	
TC Lister Scissor	MPX1980		1 pc	Can be sharpened
Curved nibbler scissor	MP1811		1 pc	
Curved nibbler scissor TC	X1811A		1 pc	Can be sharpened






Woodcast Flatbed Heater	REF	Size	Batch	Colors
Large flatbed heater for Woodcast	A001	10cm x 44cm x 88cm	1 pc	White





Keeps Woodcast sheets ready to use

Woodcast Express Heater	REF	Size	Batch	Colors
Faster heating device for Woodcast. Note: 80cm long sheets don't fit in this device unless you use the Folded variants.	A002	65cm x 23cm x 39cm	1 pc	White



Woodcast TCC-kit	REF	Size	Batch	Colors
One kit that contains everything needed for one Total Contact cast. Comes in two thicknesses. 	TCC6020	6mm thickness	1 kit	
	TCC8020	6mm thickness	1 kit	

OneStep Premade	REF	Size	Batch	Colors
Pre-molded ankle and elbow supports for paramedic use	C1465PM (Elbow)	2mm x 14.5cm x 65 cm	5 pcs	
	B1565PM (Ankle)	4mm x 15cm x 65cm	5 pcs	

Ankle brace	REF	Size	Batch	Colors
Moldable ankle brace, one size fits all. 	YL2808	One size fits all	10 pcs	Black