


Robust

Aurelia® Robust™ Blue Nitrile
Powder Free Examination Gloves



Textured for Grip
Maximum Protection

- 100% latex free, blue color
- 4.5 mil thickness for improved puncture resistance
- Micro-textured finish for improved gripping action



AURELIA®
Sign of Comfort™

www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

Robust

Product Description

Aurelia® Robust™ Blue Nitrile (Non-Latex) Powder Free Examination Gloves are made from 100% Nitrile (Non-Latex), ambidextrous, blue color, fully micro-textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com



Aurelia® Robust™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Small	26.1	550
Medium	26.5	550
Large	26.1	595
Extra Large	26.2	580
ASTM D6319-00	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Small	24.0	550
Medium	24.0	555
Large	24.0	585
Extra Large	23.7	570
ASTM D6319-00	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Small Aurelia Robust ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	4.5 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Medium Aurelia Robust ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	4.5 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Large Aurelia Robust ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	4.5 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Extra Large Aurelia Robust ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	4.5 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

ITEM NUMBER

Small	Medium	Large	Extra Large
93896	93897	93898	93899

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor: