

NEW AND IMPROVED FLUIDSHIELD* MASKS WITH SO SOFT* LINING


Where Softness And Protection Come Face-to-Face

SAME FLUIDSHIELD* PROTECTION, NEW SOFTER FEEL.


Now the number one brand of fluid-resistant face masks are softer than ever before. New and improved HALYARD* FLUIDSHIELD* Face Masks:

- Now feature our super soft SO SOFT* inner lining for more comfort than ever
- Found to be the softest and silkiest among leading surgical masks tested¹
- Powered by patented POWERGUARD* filter technology for superior microbial barrier protection
- Meets the ASTM F2100-11 Level 3 standard for the highest level of fluid-resistant protection

Four Layers Of Protection



The only masks with POWERGUARD* Filter Media for protection and SO SOFT* inner lining for comfort.



¹Halyard* Sensory Panel Study, 2013

NEW AND IMPROVED PACKAGING DESIGN



Not only has Halyard set a new performance standard with FLUIDSHIELD* Face Masks, we've created a new packaging design that makes it easy to identify product features and match the right level of protection to the task.

Features include:

- Highly visible color scheme
- Large, easy-to-read type that clearly identifies the mask's ASTM F2100-11 Level 3 protection information
- Icons communicate important information at a glance
 - Level 3 protection
 - Style of mask and color
 - Additional product features such as WrapAround visor



KNOWLEDGE NETWORK* Accredited Education
Ongoing Customer Support
Expert Sales Force
Tools & Best Practices
Clinical Research
Commitment To Excellence

For more information, please call your sales representative, or visit:

www.halyardheath.com

Call 1-844-HALYARD (1-844-425-9273) in the United States and Canada.

NOW PRODUCT SELECTION IS EASIER THAN EVER.

PROCEDURE MASKS WITH EARLOOPS

Code	Description	Items/Box	Boxes/Case
47107	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining/Orange	40	10
47137	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining, Visor/Orange	25	4
47147	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining, WrapAround Visor/Orange	25	4
47149	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining, WrapAround Visor/Orange (for use with PPE Dispensing System)	25	4
00146	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining, WrapAround Visor/Blue	25	4
00148	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining/Blue	40	10
47297	FLUIDSHIELD* Level 3 Teddy Bear Procedure Mask with SO SOFT* Lining/Teddy Bear Print	40	10
47298	FLUIDSHIELD* Level 3 Teddy Bear Procedure Mask with SO SOFT* Lining, WrapAround Visor/Teddy Bear Print	25	4

SURGICAL MASKS WITH TIES

Code	Description	Items/Box	Boxes/Case
48207	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with SO SOFT* Lining/Orange	50	6
48208	FLUIDSHIELD* Level 3 Surgical Mask with SO SOFT* Lining/Orange	50	6
48237	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with SO SOFT* Lining/Visor/Orange	25	4
48247	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with SO SOFT* Lining, WrapAround Anti-Glare Visor/Orange	25	4
48248	FLUIDSHIELD* Level 3 Surgical Mask with SO SOFT* Lining/WrapAround Visor/Orange	25	4
48297	FLUIDSHIELD* Level 3 Teddy Bear Surgical Mask with SO SOFT* Lining/Teddy Bear Print	50	6

Test Method	Test Description	LEVEL 3	
		ASTM F2100-11	FLUIDSHIELD* Masks ¹
ASTM F1862	Fluid Resistance with synthetic blood, in mm Hg	160 mm Hg	160 mm Hg
MIL-M-36954C	Delta P Differential Pressure, mm H ₂ O/cm ²	<5.0 mm H ₂ O/cm ²	3.0 mm H ₂ O/cm ²
ASTM F2101	Bacterial Filtration Efficiency (BFE), %	≥ 98%	99.8%
ASTM F2299	Submicron Particle Filtration (PFE) @ 0.1µ, %	≥ 98%	98.9%
16 CFR Part 1610	Flammability	Class 1	Class 1

¹Based on the calculated median.

