

SAFE BARIATRIC PATIENT HANDLING

Minimum Crew.

- No waiting for additional personnel or equipment when every second counts
- Increase completed calls per shift and generate more revenue due to faster transport times



2011
EMSWORLD
INNOVATION
AWARD

Maximum Safety.

- Lift, transfer, and transport bariatric patients without risking the safety of your crew or patient
- Minimize manual lifting and decrease costly workers' compensation injuries



Evacuation HoverJack® II Device

The Evacuation HoverJack® II Device is HoverTech's all-in-one transport, lift, and stair evacuation solution that minimizes many of the challenges EMS crews face in taking bariatric patients from their home to the hospital.

Evacuation HoverJack® Device Benefits:

- Air chambers inflate to lift the bariatric patient from the floor in a supine or seated position, creating a no-lift environment for the crew
- Tapered head-end, adjustable inflation chambers, and steering straps improve maneuverability around corners and down most stairwells
- Quick-release valves expedite deflation to lower patient when needed
- Support skirt with replaceable Teflon coated bottom allows device to be used across rough terrain
- Weight capacity of 700 lbs. for stairwell evacuations or 1200 lbs. for vertical lifts accommodates most patients
- Dimensions: 32" W x 78" L – HJ32EV-2

HoverMatt® Air Transfer System

The HoverMatt® Air Transfer System allows EMTs to laterally transfer patients from their bed to the Evacuation HoverJack II or to a stretcher with ease. When inflated, a cushion of air beneath the mattress reduces the force needed to move a patient by 80–90%. This reduction in force minimizes crew and patient injuries that are often caused by lateral transfers.

HoverMatt® Air Transfer System Benefits:

- Air transfer technology allows first responders to laterally transfer bariatric patients quickly and safely, while reducing patient skin shear and bruising
- Weight limit of 1200 lbs. accommodates most patients
- Available in four sizes: 28" W x 78" L – HM28HS; 34" W x 78" L – HM34HS; 39" W x 78" L – HM39HS; 50" W x 78" L – HM50HS

Double-coated and Single-Patient Use models are also available.



For more information, contact your local Horizon representative or reach us at:

877.721.5510

www.HorizonMedicalProducts.com

info@horizonmp.com



513 S. Clewell Street
Bethlehem, PA 18015
Phone 800.471.2776
Fax 610.694.9601
www.HoverMatt.com