

NEW

microdot[®]

Powder Free Nitrile Examination Gloves

SUPERIOR TOUCH SENSITIVITY AND STRENGTH



- **LATEX FREE FORMULATION**
no natural rubber protein.
- **SUPER SOFT**
and fatigue free.
- **NON – STERILE**
Single Use.
- **AMBIDEXTROUS**
- **POWDER FREE,**
no cornstarch.
- **MICRO (FINGER) TEXTURED**
surface for secure grip.
- **TESTED FOR USE**
with Chemotherapy Drugs
& Fentanyl Citrate injection.

Manufactured in compliance with the following international standards:

- ASTM D6319** Standard Specification for Nitrile Examination Gloves for Medical Application.
- ASTM F1671** Standard Test Method for Resistance of Materials used In Protection Clothing to Penetration by Blood-Borne Pathogens
- ASTM F1671** Standard for handling Chemotherapy drugs.

Size	Item No.	Packaging	
S	1201-10	100 Gloves/Box	10 Boxes/Case
M	1202-10	100 Gloves/Box	10 Boxes/Case
L	1203-10	100 Gloves/Box	10 Boxes/Case
XL	1204-10	100 Gloves/Box	10 Boxes/Case

Sold by the Case

Ask for *microdot*[®]

Ordering information

toll free: 877.374.4062 | info@microdotcs.com

NEW

microdot[®]

Powder Free Nitrile Examination Gloves

Product Specifications & Compliance:

Glove Thickness	mm	mil
Cuff	.05	2.0
Palm	.08	3.15
Finger	.09	3.54

Tensile Strength (MPa) & Elongation	Before Aging	After Aging
Tensile Strength Min. (MPa)	15	14
ASTM Requirement Min. (MPa)	14	14
Elongation Minimum Min. (%)	500	400
ASTM Requirement Min (%)	500	400

Glove Dimensions by Size	Palm Width (mm)	Total Length (mm)
S	80	230
M	95	230
L	110	230
XL	120	230



Pack

- 100 gloves per Box
- 10 boxes per Case
- 105 cases per Pallet

Ask for *microdot*[®]

Ordering information

toll free: 877.374.4062 | info@microdotcs.com

CS CAMBRIDGE SENSORS USA, LLC