

TRANSMATT

DECEDENT TRANSFER SYSTEM

Perform decedent handling tasks more efficiently while reducing risk of injury to personnel.

The TransMatt™ Decedent Transfer System allows funeral home staff to laterally transfer the deceased from their bed to a stretcher with ease. When inflated, a cushion of air beneath the mattress reduces the force needed to move a decedent by 80-90%. This reduction in force facilitates decedent handling with fewer resources, while minimizing the risk of musculoskeletal injuries to staff and promoting career longevity.



TRANSMATT™ SYSTEM BENEFITS

- Air-assisted technology:
 - Helps reduce staff injuries and associated workers' compensation costs.
 - Facilitates handling of decedents in a dignified manner and accommodates up to 1200 lbs.
 - Minimizes manual handling, while utilizing fewer staff to perform retrievals.
- Convenient, disposable product eliminates contamination concerns.
- Portability of system enables use practically anywhere.

COMPLY WITH INDUSTRY SAFETY STANDARDS

The National Institute for Occupational Safety and Health (NIOSH) recommends that no worker should lift anything heavier than 35 lbs. without using some type of assistive lift device.

PRODUCT SPECIFICATIONS

MATERIAL SPECS

Nylon fabric

Latex free

DIMENSIONS

34" W x 78" L

PRODUCT NUMBER

TM34DIS

TM34DIS PER PACKAGE

5 per box

COMPATIBLE WITH HOVERTECH AIR SUPPLIES HTAIR1200 (Variable Speed), AIR200G (800 W),





513 S. Clewell Street Bethlehem, PA 18015

Phone 800.471.2776 Fax 610.694.9601

www.HoverMatt.com

Bro TM 1015

For more information, contact HoverTech at: 800.471.2776 www.HoverMatt.com info@hovermatt.com