

**NEW**

# *microdot*<sup>®</sup> **BLK**

## **Powder Free Nitrile Examination Gloves**

SUPERIOR TOUCH SENSITIVITY AND STRENGTH

- **TESTED FOR USE**  
with Chemotherapy Drugs  
& Fentanyl Citrate injection.
- **SUPER SOFT**  
and fatigue free.
- **POWDER FREE,**  
no cornstarch.



- **LATEX FREE FORMULATION**  
no natural rubber protein.
- **MICRO (FINGER) TEXTURED**  
surface for secure grip.
- **NON – STERILE**  
Single Use.
- **AMBIDEXTROUS**

### **Manufactured in compliance with the following international standards:**

- ASTM D6319** Standard Specification for Nitrile Examination Gloves for Medical Application.
- ASTM F1671** Standard Test Method for Resistance of Materials used In Protection Clothing to Penetration by Blood-Borne Pathogens
- ASTM D6978** Standard for handling Chemotherapy drugs.

Size	Item No.	Packaging	
S	1301-10	100 Gloves/Box	10 Boxes/Case
M	1302-10	100 Gloves/Box	10 Boxes/Case
L	1303-10	100 Gloves/Box	10 Boxes/Case
XL	1304-10	100 Gloves/Box	10 Boxes/Case
XXL	1305-10	100 Gloves/Box	10 Boxes/Case

*microdot*<sup>®</sup>

Sold by the Case

**NEW**

# *microdot*<sup>®</sup> **BLK**

## Powder Free Nitrile Examination Gloves

### Product Specifications & Compliance:

Glove Thickness	mm	mil
Cuff	.05	2.0
Palm	.08	3.15
Finger	.09	3.54

Tensile Strength (MPa) & Elongation	Before Aging	After Aging
Tensile Strength Min. (MPa)	15	14
ASTM Requirement Min. (MPa)	14	14
Elongation Minimum Min. (%)	500	400
ASTM Requirement Min (%)	500	400

Glove Dimensions by Size	Palm Width (mm)	Total Length (mm)
S	80	230
M	95	230
L	110	230
XL	120	230
XXL	130	230



#### Pack

- 100 gloves per Box
- 10 boxes per Case
- 105 cases per Pallet

*microdot*<sup>®</sup>

**Horizon Medical**  
[www.horizonmedicalproducts.com](http://www.horizonmedicalproducts.com)