

microdot[®]

Infrared Forehead Thermometer

SIZE:110x51x168mm

microdot[®]

Infrared Forehead Thermometer

microdot[®]

Infrared Forehead Thermometer

QUICK START GUIDE

Make sure forehead is clean and not covered with hair, sweat, or dirt



Scan forehead at a distance of 1-5 cm



Squeeze trigger to activate and read display

microdot[®]

Infrared Forehead Thermometer

microdot[®]

Infrared Forehead Thermometer

HA-650



1 Sec 1 Second Measurement

°C/°F °C/°F Switchable

50 Mem 50 Memory Readings

Lightbulb Red Backlight Fever Alert

Waves Infrared Technology

HA-650

1 Sec 1 Second Measurement

°C/°F °C/°F Switchable

50 Mem 50 Memory Readings

Lightbulb Red Backlight Fever Alert

Waves Infrared Technology



Measurement Range:	93.2°F - 109.4°F (34.0 - 43.0°C)
Measurement Distance:	1 - 5 cm
Resolution:	0.1 °C/°F
Accuracy:	(Laboratory) ±0.4°F (95.0 - 107.6°F) / ±0.2°C (35.0 - 42.0°C) ±0.5°F (93.2 - 94.8°F) / ±0.3°C (34.0 - 34.9°C) ±0.5°F (107.8 - 109.4°F) / ±0.3°C (42.1 - 43.0°C)
Dimension:	136 x 86 x 39 mm
Operating Condition:	Temperature: 50.0 - 104.0°F (10 - 40°C) Humidity: ≤85%RH
Storage Condition:	Temperature: -4 - 131.0°F (-20 - 55°C) Humidity: ≤93%RH
Battery:	2X 1.5V AAA Batteries
Backlight:	3 color backlight

CS CAMBRIDGE SENSORS

Distributed in the USA by:
Cambridge Sensors USA, LLC
23866 W. Industrial Drive North
Plainfield, IL. 60585
toll free: 877.374.4062 • www.microdotcs.com



REF CAT No. CS-1100-01



CS_Rev.0 7/14/20



Xiamen Ants Bro Technology Co.,Ltd.
West of 4th Floor, No 5 building, Technology business establishing center, 289 Wengjiao Road, Haicang, Xiamen, China