

S.P.A. Collection

Safe. Personal. Assistance.



## Etac Flex grab rail

Flexibility and freedom

 etac®

# Etac Flex grab rail

*Innovative design that creates possibilities*

## Flexible modular grab rails

Etac Flex grab rail is the perfect choice for those who need a little extra support in the bathroom. Thanks to its flexible modular function, it is easy to adjust the handle to suit both you and your bathroom environment. The pure Scandinavian design blends in with other products and you can choose between two colours - white and grey.

## Support for good and bad days

The material of the grab rail and its well-designed shape and surface provide a stable grip even if your hands or the handles are wet. The triangular shape and the thumb support also make you feel safer.



*Etac Flex fits perfectly in the bathroom environment.*

### Etac Flex screw mount

|                              |          |
|------------------------------|----------|
| - Extension Flex 30 cm white | 81706110 |
| - Grab rail Flex 30 cm white | 81707110 |
| - Grab rail Flex 60 cm white | 81707120 |
| - Grab rail Flex 90 cm white | 81707130 |
| - Extension Flex 30 cm grey  | 81706310 |
| - Grab rail Flex 30 cm grey  | 81707310 |
| - Grab rail Flex 60 cm grey  | 81707320 |
| - Grab rail Flex 90 cm grey  | 81707330 |

### Etac Flex adhesive mount

|                              |          |
|------------------------------|----------|
| - Extension Flex 30 cm white | 81706210 |
| - Grab rail Flex 30 cm white | 81707210 |
| - Grab rail Flex 60 cm white | 81707220 |
| - Grab rail Flex 90 cm white | 81707230 |
| - Extension Flex 30 cm grey  | 81706410 |
| - Grab rail Flex 30 cm grey  | 81707410 |
| - Grab rail Flex 60 cm grey  | 81707420 |
| - Grab rail Flex 90 cm grey  | 81707430 |

**Weight** 0.35 kg per unit (0.77 lbs)

**Maximum user weight** 100 kg

**Material** Polypropylene

**Colour** White and Volcano grey



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F if so required.

*Note; To add one or more grab rail to a 30 cm, 60 cm or 90 cm combination, it is necessary to use the Etac Flex Extensions.*

## Features



### Easy to adjust

The modular system can easily be assembled into a length and shape that suits you and your bathroom.




### Secure grip

The well-designed shape and surface of the grab rail provide a stable grip even if your hands are wet.




### Creating Possibilities

You can connect as many grab rails as you want and choose the angles that suit you. This creates possibilities.



### Attach in two ways

Attach with screws or adhesive. The latter is perfect for those who cannot or do not want to drill in the bathroom.



*Etac Flex offers a wide range of positioning options when connecting grab rails. Choose the angle between  $-90^\circ$  to  $+90^\circ$  that suits your need.*





## Safe. Personal. Assistance.

*Etac S.P.A. Collection comprises Rex toilet arm support, Flex grab rail and Relax shower seat. All three have a unique function but work together to ensure your bathroom is as comfortable, attractive and safe as possible. And to ensure you can continue being yourself.*

## Lean on us. Stay yourself.

Visit **Etac.com** for more information or get in touch and we will tell you more.