



GENOA UNIVERSITY
Medicine & Surgery Faculty
Preventive Dentistry dept.
Chair: Prof. Ugo Covani

Experimentation in plaque control in the interproximal spaces using dental floss

I.D. Annamaria Genovesi
Dr. Antonio Barone
I.D. Chiara Lorenzi
Prof. Ugo Covani

January 2004

TRIAL REPORT

INTRODUCTION

Complete domestic dental hygiene is guaranteed by daily cleaning of all dental surfaces: the treatment, thus, regards both the vestibular and lingual surfaces and the interproximal surfaces. For the plane surfaces of the teeth, one uses a toothbrush, which cannot reach the interdental zones, unless they are extremely distant. Hygiene in these zones is guaranteed by the use of specific auxiliary, which enables the removal of food residue and plaque: dental floss.

TRIAL HYPOTHESIS

Riser Floss®, *there's only one!* dental floss, whose section increases during use and, specifically, when subjected to mechanical solicitation, delivers better interdental plaque removal and helps reduce gingival bleeding.

TRIAL OBJECTIVES

Verification of the efficacy of **Riser Floss®**, *there's only one!* dental floss by comparison with other marketed dental floss types.

TRIAL MATERIALS & METHODS

The supposed efficacy of **Riser Floss®**, *there's only one!* TRIAL was verified by comparison with three of the widest used and marketed types of dental floss, in order to verify which was best suited to the needs of the patients examined.

The trial involved 40 patients divided into 4 groups.

The first group used a waxed dental floss (green container), with the following technical characteristics: *Regular*, PA 6.6, Dtex 940/136, similar to Oral B traditional waxed dental floss.

The second group used a monofilament dental floss (light blue container), with the following technical characteristics: PTFE, Dtex 1200, folded, slightly waxed similar to Gore Glide Floss.

The third used a fine dental floss (neutral container), with the following technical characteristics: *Fine*, PA 6.6, Dtex 700/104), similar to Butler Fine unwaxed floss.

Finally, the fourth group used **Riser Floss®**, *there's only one!* floss (white container), with the following technical characteristics: PL, Dtex 1300/376.

The following parodontal indices were measured in each patient at the beginning and end of the trial:

- ✚ the bacterial plaque index, applying O'LEARY's (PCR) technique, using a bi-tonal plaque indicator;
- ✚ the bleeding index, applying the (GBI) AIMANO and BAY technique.

The patients had not professional oral hygiene care but were instructed in oral hygiene technique with the use of toothbrush and dental floss.

The patients were re-examined after two weeks.

All patients had integral anatomical papillary structure. They were not informed of the type of dental floss used.

TRIAL RESULTS

All 40 patients completed the trial. The full details of the plaque and bleeding indices are provided in the table below.

For the **green container** (*Regular* waxed dental floss), at the first examination, an average PCR value of 41.60% and an average GBI value of 31.80% were recorded; while at the second examination the average PCR value was 19.90% and the average GBI value was 15.97%.

The mean percentage reduction in the plaque index was 52.16%, while the bleeding index fell by 49.78%.

For the **light blue container** (PTFE monofilament dental floss), at the first examination, an average PCR value of 53.17% and an average GBI value of 39.03% were recorded; while at the second examination the average PCR value was 29.79% and the average GBI value was 29.74%.

The mean percentage reduction in the plaque index was 43.97% while the bleeding index fell by 25.08%.

For the **neutral container** (*Fine* unwaxed dental floss) at the first examination, an average PCR value of 38.09% and an average GBI value of 24.66% were recorded; while at the second examination the average PCR value was 19.24% and the average GBI value was 16.82%.

The mean percentage reduction in the plaque index was 49.49 while the bleeding index fell by 31.79%.

For the **white container** (*Riser Floss®*, *there's only one!* dental floss) at the first examination, an average PCR value of 53.49% and an average GBI value of 18.49%; were recorded; while at the second examination the average PCR value was 16.12% and the average GBI value was 9.13%.

The mean percentage reduction in the plaque index was 69.86% while the bleeding index fell by 50.62%.

TRIAL CONCLUSIONS

All patients showed clinical improvement from the first to the second examination two weeks later, independent of the type of dental floss used.

The inter-group differences are significant, in particular with reference to the plaque index: the use of *Riser Floss®*, *there's only one!* dental floss delivered a mean reduction of almost 70%, notably superior to other flosses.

The reduction in the bleeding index was minimal only with reference to the *Regular* dental floss 50.62% and 49.78%, respectively.

Further, the *Riser Floss®*, *there's only one!* dental floss was more efficacious with respect to the other dental flosses in the removal of dental plaque in patients with moderate to severe parodontitis.

TABLES

Neutral container (Fine unwaxed dental floss)

| Pos. | Age | Sex | PCR1 | PCR2 | GBI1 | GBI2 |
|----------------|------------|------------|---------------|---------------|---------------|---------------|
| 01 | 55 | F | 35,52% | 15,78% | 17,10% | 5,26% |
| 02 | 56 | F | 60,71% | 46,43% | 50,00% | 46,43% |
| 03 | 42 | M | 41,34% | 20,19% | 29,80% | 25,00% |
| 04 | 22 | M | 17,19% | 4,69% | 9,37% | 3,12% |
| 05 | 21 | F | 66,07% | 16,07% | 24,10% | 7,14% |
| 06 | 47 | F | 44,23% | 33,65% | 33,65 | 16,35% |
| 07 | 47 | M | 56,25% | 35,71% | 40,18% | 32,14% |
| 08 | 29 | F | 16,38% | 12,93% | 19,83% | 1,78% |
| 09 | 21 | M | 38,39% | 5,36% | 8,03% | 1,78% |
| 10 | 26 | F | 4,84% | 1,61% | 14,52% | 13,71% |
| Average | | | 38,09% | 19,24% | 24,66% | 16,82% |

Light blue container (slightly waxed PTFE monofilament dental floss)

| Pos. | Age | Sex | PCR1 | PCR2 | GBI1 | GBI2 |
|----------------|------------|------------|---------------|---------------|---------------|---------------|
| 01 | 72 | M | 66,18% | 45,59% | 54,41% | 47,06% |
| 02 | 34 | F | 49,07% | 35,18% | 30,55% | 21,30% |
| 03 | 41 | F | 46,87% | 37,50% | 30,21% | 22,92% |
| 04 | 28 | F | 32,14% | 25,00% | 22,32% | 16,96% |
| 05 | 46 | M | 80,36% | 40,18% | 32,14% | 0,00% |
| 06 | 48 | F | 55,77% | 22,11% | 46,15% | 37,50% |
| 07 | 51 | M | 56,00% | 45,00% | 45,00% | 42,00% |
| 08 | 60 | F | 100,00% | 27,42% | 100,00% | 87,90% |
| 09 | 28 | F | 20,31% | 9,37% | 9,37% | 2,34% |
| 10 | 48 | F | 25,00% | 10,58% | 20,19% | 14,42% |
| Average | | | 53,17% | 29,79% | 39,03% | 29,24% |

Green container (Regular waxed dental floss)

| Pos. | Age | Sex | PCR1 | PCR2 | GBI1 | GBI2 |
|----------------|------------|------------|---------------|---------------|---------------|---------------|
| 01 | 68 | M | 67,19% | 45,31% | 53,12% | 41,19% |
| 02 | 27 | M | 20,37% | 12,96% | 11,11% | 4,63% |
| 03 | 25 | F | 17,50% | 9,17% | 4,17% | 0,00% |
| 04 | 32 | F | 15,52% | 10,34% | 6,90% | 4,31% |
| 05 | 67 | F | 52,94% | 39,70% | 61,76% | 52,94% |
| 06 | 52 | M | 68,52% | 13,89% | 32,41% | 8,33% |
| 07 | 20 | M | 13,71% | 9,68% | 8,87% | 6,45% |
| 08 | 21 | F | 46,09% | 9,37% | 30,47% | 5,47% |
| 09 | 64 | F | 80,88% | 29,41% | 80,88% | 20,59% |
| 10 | 25 | M | 33,33% | 19,17% | 28,33% | 15,83% |
| Average | | | 41,60% | 19,90% | 31,80% | 15,97% |

White container (*Riser Floss®*, there's only one!)

| Pos. | Age | Sex | PCR1 | PCR2 | GBI1 | GBI2 |
|----------------|-----|-----|---------------|---------------|---------------|--------------|
| 01 | 29 | F | 28,12% | 17,19% | 20,31% | 14,84% |
| 02 | 46 | M | 54,81% | 17,31% | 49,04% | 10,58% |
| 03 | 53 | M | 72,22% | 52,77% | 44,44% | 27,77% |
| 04 | 45 | M | 18,75% | 0,00% | 7,14% | 3,57% |
| 05 | 24 | F | 9,48% | 1,71% | 11,21% | 2,59% |
| 06 | 34 | F | 36,16% | 3,91% | 21,87% | 6,25% |
| 07 | 32 | F | 13,33% | 7,50% | 4,16% | 3,33% |
| 08 | 48 | F | 100,00% | 39,13% | 19,56% | 15,22% |
| 09 | 36 | F | 9,37% | 2,34% | 0,00% | 0,00% |
| 10 | 31 | F | 33,93% | 19,64% | 7,14% | 7,14% |
| Average | | | 53,49% | 16,12% | 18,49% | 9,13% |

Summarizing table

| Pos. | Type of floss | Container | PCR1 | PCR2 | GBI1 | GBI2 |
|------|---------------------------|------------|--------|--------|--------|--------|
| 01 | <i>Fine</i> unwaxed | neutral | 38,09% | 19,24% | 24,66% | 16,82% |
| 02 | <i>Ptfe</i> lightly waxed | light blue | 53,17% | 29,79% | 39,03% | 29,24% |
| 03 | <i>Regular</i> waxed | green | 41,60% | 19,90% | 31,80% | 15,97% |
| 04 | Riser floss® | white | 53,49 | 16,12% | 18,49% | 9,13% |

Average reduction in percentage indexes table

| Pos. | Type of floss | Container | PCR | GBI |
|------|---------------------------|------------|--------|--------|
| 01 | <i>Fine</i> unwaxed | neutral | 49,49% | 31,79% |
| 02 | <i>Ptfe</i> lightly waxed | light blue | 43,97% | 25,08% |
| 03 | <i>Regular</i> waxed | green | 52,16% | 49,78% |
| 04 | Riser floss® | white | 69,86 | 50,62% |

Observations

The mean result obtained with *Riser Floss®*, *there's only one!* dental floss in plaque reduction is superior by 43.92% to the average of other dental flosses.

The mean result obtained with *Riser Floss®*, *there's only one!* dental floss in the reduction of the bleeding index is superior by 42.39% to the average of other dental flosses.