



GluStitch Inc.



GluStitch Inc.

#307-7188 Progress Way  
Delta, BC Canada V4G1M6  
info@glustitch.com  
1-877-940-2262

# Designed For Dermatologists And Plastic Surgeons

GluStitch Inc Family of Products



Like us on Facebook at  
[www.facebook.com/dermglu](http://www.facebook.com/dermglu)



Shop online at  
[www.dermglu.com](http://www.dermglu.com)

## Testimonial and Catalog

# Testimonial



My name is **Dr. Shawn Allen** and I am a board certified dermatologist and fellowship trained Mohs and cosmetic surgeon in Boulder, Colorado. I am happy to talk about my experience with the GluStitch® product line. I have had great success with these products since 2009; using them in wide range of procedures from face-lifts and standard excisions to closing Mohs wounds whether they are being closed by primary or secondary intent.

I have found these products to really change my practice in some significant ways; with respect to patient care and convenience as well as improving the efficiency of my practice. Using GluSeal® and GluStitch® Twist tissue adhesives I am able

to close wounds that would normally require epidermal sutures and the substitution of these sutures with tissue adhesives allows for a quicker closure as well as a cleaner closure without the risk of suture track marks. I believe that as these adhesives polymerize along the incision line, they work to evert the wound edges which yields ideal cosmetic results. The use of adhesives versus sutures seals the wound so the risk of infection likely is decreased. Another benefit is that the patients don't need to follow a complicated wound care regimen as it is not required to apply any ointments to the surgical site or change dressings. It also has been a cost effective approach for my practice, not only do we save the time not needing to schedule the suture removal visits but the cost of GluStitch® products compare favorably to placing top sutures or using some of the other cyanoacrylate products on the market.

## Shawn Allen, MD

Director and Founder of Dermatology Specialists of Boulder  
2880 Folsom Steet Suite 200 Boulder, CO 80304

# Before and After Photos




CLOSURE OF LEFT CHEEK



18 MONTHS LATER



CLOSURE ON CHEST



6 YEARS LATER

# The GluStitch® Family of Products Catalog

GluStitch Inc manufacturers all of its products using butyl cyanoacrylate and blends of butyl and octyl cyanoacrylate. Our most popular product is the GluSeal High Viscosity 5 mL Clear Kit because practitioners can reduce their costs to as little as \$3 per application. All products are FDA cleared and come with a 100% satisfaction guarantee.



## GluSeal® High Viscosity Clear with 50 disposable pipettes and autoclavable tray

GluSeal® High Viscosity is a blend of butyl and octyl cyanoacrylate that is effective as a liquid dressing to cover and support closure of minor lacerations. The combination of butyl and octyl forms a strong yet flexible dressing to protect wounds. The thicker formula improves the ease of application.



## GluSeal® Single-Use

GluSeal® Single-use is available in clear or violet formulations. This presentation comes with 12 x 0.2 mL clip-and-use prefilled ampoules.



## GluStitch® Twist

GluStitch® Twist is made available in 6 packs of 0.5mL ampoules that pouched individually. These sterile units can be purchased through our exclusive distributor SkinStitch Corp based in Massena, NY. 1-866-254-3900

Contact dealer for pricing information