

## Bion, Comfy Pushchair

BION Code: 00643



<b>Seat Width</b>	<b>16½" (42cm)</b>
<b>Seat Depth</b>	<b>16" (41cm)</b>
<b>Seat Height</b>	<b>17¾" (45cm)</b>
<b>Backrest Height (From Seat)</b>	<b>15¼" (39cm)</b>
<b>Distance between Armrest</b>	<b>17¾" (45cm)</b>
<b>Overall Width</b>	<b>23" (58cm)</b>
<b>Overall Height</b>	<b>33¾" (86cm)</b>
<b>Overall Length</b>	<b>40½" (103cm)</b>
<b>Net Weight</b>	<b>12kg</b>
<b>Weight Capacity</b>	<b>100kg</b>
<b>Front Wheel Size</b>	<b>6"</b>
<b>Rear Wheel Size</b>	<b>12" Solid Tyres, Mag Rim</b>
<b>Frame Material</b>	<b>Aluminium</b>
<b>Armrest</b>	<b>Half-Length, PU, Fixed</b>
<b>Footrest</b>	<b>Composite Footplate, Flip-Up</b>
<b>Foldable Backrest</b>	<b>✓</b>
<b>Removable Backrest Cover</b>	<b>✓</b>
<b>Removable Seat Cover</b>	<b>✓</b>
<b>User Brakes</b>	<b>✓</b>
<b>Attendant Brakes</b>	<b>✓</b>
<b>Safety Belt</b>	<b>✓</b>
<b>Anti-Tipper</b>	<b>Optional</b>

*1 year warranty on frame only. Wear and tear parts are not covered under warranty.*