

Bion, Comfy Wheelchair

BION Code: 00644



Seat Width	16½" (42cm)
Seat Depth	16" (41cm)
Seat Height	18½" (47cm)
Backrest Height (From Seat)	15" (38cm)
Distance between Armrest	17¾" (45cm)
Overall Width	25" (64cm)
Overall Height	35" (89cm)
Overall Length	40½" (103cm)
Net Weight	14kg
Weight Capacity	100kg
Front Wheel Size	7"
Rear Wheel Size	22" x 1⅜" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Half-Length, PU, Fixed
Footrest	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.