



## Bion, Rollator, Aluminium, w Seat

BION Code: 00834



<b>Overall Height</b>	<b>78cm - 89cm</b>
<b>Inner Width between Handles</b>	<b>45cm</b>
<b>Seat Width</b>	<b>36cm</b>
<b>Seat Depth</b>	<b>33cm</b>
<b>Seat Height</b>	<b>53cm</b>
<b>Backrest Height</b>	<b>26cm</b>
<b>Weight</b>	<b>7.6kg</b>
<b>Weight Capacity</b>	<b>90kg</b>
<b>Frame Material</b>	<b>Aluminium</b>

*1 year warranty on frame only. Wear and tear parts are not covered under warranty.*

*Updated as of 11/12/2021*