

Bion, iLight Pushchair, Detachable, 18"

BION Code: 01224



Seat Width	18" (46cm)
Seat Depth	17" (43cm)
Seat Height	18" (46cm)
Backrest Height (From Seat)	16¼" (42cm)
Distance between Armrest	18½" (47cm)
Overall Width	24" (61cm)
Overall Height	33½" (85cm)
Overall Length	40½" (103cm)
Net Weight	13kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	12" Solid Tyres, Mag Rim
Frame Material	Aluminium
Armrest	Half-Length, Padded, Flip-Up
Footrest	Composite Footplate, Detachable
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.