


bion

Bion, Anti-Bedsore Mattress, V250

Item code: 02002




prevention & treatment of mild pressure sores*



alternating inflation method



adjustable firmness



low vibration



low noise



waterproof mattress





Motor Pump Specifications

Size (L x W x H):	27 x 12.5 x 8 cm
Weight:	1.35 kg
Material:	ABS Plastic
Power Supply:	220-240V, 50Hz 15VA
Fuse Rating:	T500mA, 250V
Cycle Time:	12 mins
Pressure Range:	30-130±10mmHg
Environment:	Atmospheric Pressure
	80KPa to 106KPa
	Temperature
	Operation: 5°C to 35°C (41°F to 95°F) Storage: -20°C to 55°C (-4°F to 131°F) Shipping: -20°C to 55°C (-4°F to 131°F)

Humidity	Operation: 15% - 85%
	Storage: 10% to 90%RH, but not requiring a water vapour partial pressure greater than 50hPa
	Shipping: 10% to 90%RH, but not requiring a water vapour partial pressure greater than 50hPa

Product Diagram



Pressure Relief Mattress Specifications

Dimension (L x W x H):	189 x 90 x 7 cm
Weight:	1.95 kg
Max. Support Weight	135 kg
Cell Material	PVC
Cover	-

Scan QR code for more information on product and warranty

