

## Bion Shower Chair 797

Item code: 02176





### Support for Standing and Sitting

Bion Shower Chair 797 is basic shower chair with 4 legged anti-slip stoppers and height adjustable legs. Side handles are available to help in standing and sitting.

### Specifications:

Seat Width	15.2" / 38.5cm
Seat Depth	11.6" / 29.5cm
Seat Height	13.8"-18.1" / 35cm-46cm
Backrest Height (From Seat)	NIL
Distance between Armrests	NIL
Overall Width	19.7" / 50cm
Overall Height	15.4"-19.3" / 39cm-49cm
Overall Length	15.7" / 40cm
Net Weight	1.75kg
Weight Capacity	90kg
Frame Material	Aluminium
Armrests	X

### Technical Drawings:



Scan QR code for more information on product and warranty