

Bion, iLight Wheelchair, EZ

BION Code: 02286



Seat Width	17 ³ / ₄ " (45cm)
Seat Depth	17" (43cm)
Seat Height	18 ¹ / ₂ " (47cm)
Backrest Height (From Seat)	15 ³ / ₄ " (40cm)
Distance between Armrest	17" (43cm)
Overall Width	26" (66cm)
Overall Height	35" (89cm)
Overall Length	40 ¹ / ₂ " (103cm)
Net Weight	12.8kg
Weight Capacity	100kg
Front Wheel Size	7"
Rear Wheel Size	22" x 1 ³ / ₈ " Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Half-Length, Padded, Fixed
Footrest	Composite Footplate, Fixed
Foldable Backrest	✓
Removable Backrest Cover	✗
Removable Seat Cover	✗
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.