

Bion, Commode, Stationary, 103



Item code: 02295



Specifications:

Overall Width	20.9" / 53cm
Overall Height	35"-39.3" / 89-100cm
Overall Length	19.7" / 50cm
Distance between Armrest to Armrest	17.1" / 43.5cm
Seat Width	15.4" / 39cm
Seat Depth	16.5" / 42cm
Seat Height	17.3"-21.3" / 44-54cm
Backrest Height	18.5" / 47cm
Clearance for Toilet Bowl	13.2"-17.3" / 33.5-44cm
Net Weight	4.4kg
Weight Capacity	100kg
Wheels Diameter	NIL
Number of Lockable Wheels	NIL
Armrests Type	ABS Armrest, Fixed
Footrests	NIL
Frame Material	Aluminium
Foldable Frame	✓
Seat Cover	✓
Toilet Bowl Friendly	✓
Rear Anti-Tipper	X
Front Anti-Tipper	X
Back Lift	X
Safety Belt	Optional

Technical Drawings



Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty



Scan QR code for more information