

bion

Bion Traveller

Item code: 02297




Lightweight & Compact

Travel around the world with Bion Traveller, your sturdy, lightweight, and compact travel companion. It comes with a convenient travel bag.



Dimensions:



Specifications:

Seat Width	14.2" / 36cm
Seat Depth	14.8" / 37.5cm
Seat Height	18.5" / 47cm
Backrest Height (From Seat)	16.1" / 41cm
Distance between Armrest	14.8" / 37.5cm
Overall Width (Open)	16.9" / 43cm
Overall Height	34.4" / 87.5cm
Overall Length (Open)	33" / 83.75cm
Net Weight	7.3kg
Weight Capacity	80kg
Front Wheel Size	5"
Rear Wheel Size	8" Solid Tyres, Mag Rim
Frame Material	Aluminium
Armrests	Half-Length, Rubber, Flip-Up
Footrests	Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	X
Removable Seat Cover	X
User Brakes	X
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	X

Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty