# bion

## Bion iLight Detachable Wheelchair 16"

Item code: 02310



## Lightweight, Full-Featured

Bion iLight Detachable Wheelchair is part of our Full-Featured Series; designed for patients. Flip-up armrests and removable leg rests allow easy patient transfer from bed to wheelchair.

Removable cushions are easy to change and machine wash. Antitippers prevent the wheelchair to tip backward when travelling up-slope.



### **Specifications:**


Seat Width Seat Depth **Seat Height Backrest Height** (From Seat) Distance between Armrest Overall Width (Open) **Overall Height** Overall Length (Open) **Net Weight Weight Capacity Front Wheel Size Rear Wheel Size** Frame Material **Armrests Footrests** Foldable Backrest


18.9" / 48cm
15.7" / 40cm
16.5" / 42cm
22.8" / 58cm
35" / 89cm
41" / 104cm
14.8kg
100kg
7"
22" X 13/8"
Solid Tyres, Spokes Rim
Aluminium
Half-Length, PU, Fixed
Composite Footplate,
Detachable

15.4" / 39cm

16.9" / 43cm

#### **Dimensions:**





#### **Warranty Coverage**

Removable Backrest Cover
Removable Seat Cover
User Brakes

Attendant Brakes
Safety Belt

**Anti-Tipper** 

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty

28.7" (73cm)