

Bion, Comfy Pushchair, 3G

BION Code: 02311



Seat Width	15¾" (40cm)
Seat Depth	15¾" (40cm)
Seat Height	17½" (44.5cm)
Backrest Height (From Seat)	16¼" (41.5cm)
Distance between Armrest	18" (45.5cm)
Overall Width	21½" (54.5cm)
Overall Height	34½" (87.5cm)
Overall Length	39" (99cm)
Net Weight	10.4kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	12" Solid Tyres, Mag Rim
Frame Material	Aluminium
Armrest	Half-Length, Padded, Fixed
Footrest	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✗
Removable Seat Cover	✗
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.