

# bion

## Bion Comfy Wheelchair 3G

Item code: 02316



### Comfortable & Compact


Bion Comfy Wheelchair 3G is part of our Compact Series. With a foldable footrest and backrest, the unit can be folded into a more portable form. In addition to being compact, this unit has PVC seat covers which makes wiping and cleaning easier.



### Specifications:

Seat Width	15.7" / 40cm
Seat Depth	15.7" / 40cm
Seat Height	17.9" / 45.5cm
Backrest Height (From Seat)	16.5" / 42cm
Distance between Armrest	17.9" / 45.5cm
Overall Width (Open)	24" / 61cm
Overall Height	35" / 89cm
Overall Length (Open)	39" / 99cm
Net Weight	11.5kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	20" X 1 3/8" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrests	Half-Length, Padded, Fixed
Footrests	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	X
Removable Seat Cover	X
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

### Dimensions:



### Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty