

Bion, Comfy Wheelchair, 3G

BION Code: 02316



Seat Width	15¾" (40cm)
Seat Depth	15¾" (40cm)
Seat Height	17¼" (45.5cm)
Backrest Height (From Seat)	16½" (42cm)
Distance between Armrest	18" (45.5cm)
Overall Width	24" (61cm)
Overall Height	35" (89cm)
Overall Length	39" (99cm)
Net Weight	11.5kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	20" x 1⅜" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Half-Length, Padded, Fixed
Footrest	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✗
Removable Seat Cover	✗
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.