

Bion, iLight Pushchair, EZ

BION Code: 02318



Seat Width	17¾" (45cm)
Seat Depth	16¾" (43cm)
Seat Height	18" (46cm)
Backrest Height (From Seat)	15¼" (39cm)
Distance between Armrest	16¾" (43cm)
Overall Width	25" (65cm)
Overall Height	33½" (85cm)
Overall Length	40½" (103cm)
Net Weight	10.65kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	12" Solid Tyres, Mag Rim
Frame Material	Aluminium
Armrest	Half-Length, Padded, Fixed
Footrest	Composite Footplate, Fixed
Foldable Backrest	✓
Removable Backrest Cover	✗
Removable Seat Cover	✗
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.