

Bion, iLight Wheelchair, Detachable, HD, 20"

BION Code: 02357



Seat Width	20" (51cm)
Seat Depth	18" (46cm)
Seat Height	20" (51cm)
Backrest Height (From Seat)	16" (41cm)
Distance between Armrest	20½" (52cm)
Overall Width	27½" (70cm)
Overall Height	36½" (93cm)
Overall Length	42¾" (109cm)
Net Weight	19kg
Weight Capacity	120kg
Front Wheel Size	8"
Rear Wheel Size	24" x 1⅜" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Half-Length, PU, Flip-Up
Footrest	Composite Footplate, Detachable
Foldable Backrest	*
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.