

Bion, iLight Wheelchair, Elevating

BION Code: 02365



Seat Width	17½" (44cm)
Seat Depth	16¾" (43cm)
Seat Height	18½" (47cm)
Backrest Height (From Seat)	15½" (39cm)
Distance between Armrest	17" (43cm)
Overall Width	26½" (67cm)
Overall Height	35" (89cm)
Overall Length - Legrest Elevated	45" (114cm)
	47" (119cm)
Net Weight	15.5kg
Weight Capacity	100kg
Front Wheel Size	7"
Rear Wheel Size	22" x 1⅜" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Full Length, Padded, Flip-Up
Footrest	Composite Footplate, Detachable, Elevating
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.