

Bion, Shower Chair, w Armrests 004

BION Code: 02383



Seat Width	16¼" (41.5cm)
Seat Depth	15¼" (38.5cm)
Seat Height	14"-18" (35.5-45.5cm)
Backrest Height (From Seat)	15" (38cm)
Distance between Armrests	18" (46cm)
Overall Width	21" (53.5cm)
Overall Height	28½"-32¾" (72.5cm-83cm)
Overall Length	22" (56cm)
Net weight	3.4kg
Weight Capacity	90kg
Frame Material	Aluminium
Armrests	Yes

1 year warranty on frame only. Wear and tear parts are not covered under warranty.